

SurePrep™ Urine Bacterial RNA Purification Kit

Product Cat. # BP2804-50

Instruction Manual

I.	Introduction	2
	A. Product Description	
	B. Overview of Procedure	
	C. Kit Specifications	
	D. Advantages	
	E. Kit Components	
	F. Storage Conditions and Product Stability	
	G. Precautions and Disclaimers	
	H. Customer-Supplied Reagents and Equipment	
	I. Working with RNA	
II.	Set-Up and Preparation of Sample Lysate	6
	A. Equipment Preparation	
	B. Preparation of Lysate from Leukocytes	
III.	Purifying Total RNA from Sample Lysate	8
	A. Binding RNA to Column	
	B. Column Wash	
	C. RNA Elution	
	D. Assessing RNA Yield and Quality by UV Absorbance	
	E. Storage of RNA	
IV.	Troubleshooting Guide	10
V.	Related Products	12
	A. Additional RNA Purification Kits	
	B. Other Fisher BioReagents Functionally Tested for RNA Research	
VI.	Appendix A	13
	Protocol for Optional On-Column DNA Removal	

I. Introduction

A. Product Description

The SurePrep™ Urine Bacterial RNA Purification Kit is designed for the rapid preparation of bacterial RNA from urine samples. Bacterial RNA can be isolated from both human urine samples and urine samples from animals in order to study the levels and types of bacteria that are present, as well as to study the stage of bacterial pathogenesis through the use of RNA biomarkers. The kit allows for the isolation of RNA from both Gram negative and Gram positive bacteria, including *E. coli*, *Proteus* spp., *Klebsiella* spp., *Enterobacter* spp., *Serratia* spp., *Pseudomonas* spp, *Clostridial* spp. and *Leptospirosis* spp., as well as *Chlamydia trachomatis* and *Neisseria gonorrhoeae*. Fisher's Urine Bacterial RNA Purification Kit purifies all sizes of RNA, from large mRNA and ribosomal RNA down to microRNA (miRNA) and small interfering RNA (siRNA). The RNA is preferentially purified from other cellular biomolecules such as proteins, as well as from the contaminating chemical species found in urine such as glucose and salts, without the use of phenol or chloroform. Typical yields of RNA will vary depending on the urine sample and the bacterial species, if any, present in the urine. The purified RNA is of the highest integrity, and can be used in a number of downstream applications including real time PCR, reverse transcription PCR, Northern blotting, RNase protection and primer extension, and expression array assays.

B. Overview of Procedure

Purification is based on spin column chromatography using a proprietary resin as the separation matrix. The RNA is preferentially purified from other cellular components such as proteins without the use of phenol or chloroform. The process involves first pelleting the bacteria that are present in the urine sample through the use of centrifugation, followed by lysing the bacterial cells with the provided Lysis Solution (please see the flow chart on page 5). Ethanol is then added to the bacterial lysate, and the solution is loaded onto a spin-column. Fisher's resin binds RNA in a manner that depends on ionic concentrations. Thus only the RNA will bind to the column, while the contaminating proteins will be removed in the flowthrough or retained on the top of the resin. The bound RNA is then washed twice with the provided Wash Solution in order to remove any remaining impurities, and the purified total RNA is eluted with the Elution Buffer.

C. Kit Specifications

Kit Specifications	
Column Binding Capacity	50 µg
Volume of Urine Processed	10 – 50 mL
Size of RNA Purified	All sizes, including small RNA (<200 nt)
Time to Complete 10 Purifications	30 minutes
Average Yield	~ 0.5 µg RNA per 1 x 10 ⁷ cells (Varies due to cell density of sample)

D. Advantages

- Fast and easy processing using rapid spin-column format
- Isolate total RNA, from large mRNA down to microRNA (miRNA)
- Isolate high quality total RNA from urine
- No phenol or chloroform extractions

E. Kit Components

Component	Catalog # BP2804-50 (50 preps)
Lysis Solution	40 mL
Wash Solution	22 mL
RNA Elution Buffer	6 mL
Spin Columns	50
Collection Tubes	50
Elution tubes (1.7 mL)	50
Product Insert	1

F. Storage Conditions and Product Stability

All solutions should be kept tightly sealed and stored at room temperature. These reagents should remain stable for at least 2 years in their unopened containers.

G. Precautions and Disclaimers

This kit is designed for research purposes only. It is not intended for human or diagnostic use.

Ensure that a suitable lab coat, disposable gloves and protective goggles are worn when working with chemicals. For more information, please consult the appropriate Material Safety Data Sheets (MSDSs). The MSDS can be requested through our Customer Service.

H. Customer-Supplied Reagents and Equipment

You must have the following in order to use the SurePrep™ Urine Bacterial RNA Purification Kit:

- Benchtop microcentrifuge
- Swinging bucket centrifuge
- 96 – 100% ethanol
- β -mercaptoethanol
- 50 mL conical tubes
- Microcentrifuge tubes
- Lysozyme-containing TE Buffer (3 mg/mL lysozyme in TE Buffer)


I. Working with RNA

RNases are very stable and robust enzymes that degrade RNA. Autoclaving solutions and glassware is not always sufficient to actively remove these enzymes. The first step when preparing to work with RNA is to create an RNase-free environment. The following precautions are recommended as your best defense against these enzymes.

- The RNA area should be located away from microbiological work stations.
- Clean, disposable gloves should be worn at all times when handling reagents, samples, pipettes, disposable tubes, etc. It is recommended that gloves are changed frequently to avoid contamination.
- There should be designated solutions, tips, tubes, lab coats, pipettes, etc. for RNA only.
- All RNA solutions should be prepared using at least 0.05% DEPC-treated autoclaved water or molecular biology grade nuclease-free water.
- Clean all surfaces with commercially available RNase decontamination solutions.
- When working with purified RNA samples, ensure that they remain on ice during downstream applications.

Flow Chart

Procedure for Purifying Bacterial RNA using Fisher's Urine Bacterial RNA Purification Kit



II. Set-Up and Preparation of Sample Lysate

A. Equipment Preparation

Centrifugation steps are carried out in both a benchtop microcentrifuge and a swinging bucket centrifuge. Various speeds are required for different steps, so please check your microcentrifuge specifications to ensure that it is capable of the proper speeds. All centrifugation steps are performed at room temperature. The correct rpm can be calculated using the formula:

$$RPM = \sqrt{\frac{RCF}{(1.118 \times 10^{-5}) (r)}}$$

where *RCF* = required gravitational acceleration (relative centrifugal force in units of g); *r* = radius of the rotor in cm; and *RPM* = the number of revolutions per minute required to achieve the necessary *g*-force.

If you do not own a variable speed microcentrifuge consider purchasing Thermo Scientific's Sorvall Legend benchtop model that offers power, safety and convenience.

- Choice of 17,000 or 21,000 x g (for RNA purification using SurePrep kits, the microcentrifuge with 17,000 x g is sufficient)
- Holds 36 x 0.5 mL microtubes, 24 x 2 mL tubes or 8 x 8 PCR
- Unique ClickSeal™ bio-containment rotor lid for safe processing of infectious specimens
- Fast acceleration and deceleration speeds up your protocols
- Broad range of rotors supports virtually any application
- Intuitive controls and vivid display
- Highly resistant materials allow vigorous cleaning and autoclaving

Sorvall Legend Micro Centrifuges

Technical Specifications

	Sorvall Legend Micro 17 & 17R	Sorvall Legend Micro 21 & 21R
Max g-force:	17,000	21,100
Max RPM:	13,300	14,800
Noise level:	<55 dBA	<56 dBA
Time set range:	1 min - 99 min; 1 min increments	1 min - 99 min; 1 min increments
Temp set range:	Set from -9 °C to +40 °C; per 1 °C increment	Set from -9 °C to +40 °C per 1 °C increment

Ordering Information

	Cat. No.	Cat. No.
Sorvall Legend Micro 17/17R	230V 50/60Hz	120V 60 Hz
Sorvall Legend Micro 17, includes 24 x 1.5/2.0 mL rotor with ClickSeal bio-containment lid	75002430	75002431
Sorvall Legend Micro 17R, includes 24 x 1.5/2.0 mL rotor with ClickSeal bio-containment lid	75002440	75002441
Sorvall Legend Micro 21/21R	230V 50/60Hz	120V 60 Hz
Sorvall Legend Micro 21, incl. 24 x 1.5/2.0 mL rotor with ClickSeal bio-containment lid	75002435	75002436
Sorvall LegendMicro 21R, incl. 24 x 1.5/2.0 mL rotor with ClickSeal bio-containment lid	75002445	75002446

For detailed product specifications, information on additional rotors, lids and adapters visit www.thermo.com

B. Preparation of Lysate from Urine Bacteria

Notes Prior to Use

- It is recommended that at least 20 – 30 mL of urine be used as the input for each column. The maximum urine input should not exceed 50 mL.
- Ensure that all solutions are at room temperature prior to use, and that no precipitates have formed. If necessary, warm the solutions and mix well until the solutions become clear again.
- Prepare a working concentration of the **Wash Solution** by adding 50 mL of 95% ethanol (provided by the user) to the supplied bottle containing the concentrated **Wash Solution**. This will give a final volume of 72 mL. The label on the bottle has a box that may be checked to indicate that the ethanol has been added.
- Prepare an appropriate amount of Lysis Solution by adding 10 μ L of β -mercaptoethanol (provided by the user) to each 1 mL of Lysis Solution required. β -mercaptoethanol is toxic and should be dispensed in a fume hood.
- Prepare an appropriate amount of TE Buffer containing 3 mg/mL of lysozyme. This solution should be prepared with sterile, RNase-free TE Buffer and RNase-free lysozyme. Keep the solution on ice until needed. These reagents are to be provided by the user.
- The use of fresh urine is recommended for this procedure. Since RNA is a biomolecule that undergoes rapid turnover within a bacterial cell, a delay in isolation may result in purified RNA that does not reflect its true state at the time of sample collection. Alternatively, the bacteria can be pelleted, flash-frozen immediately using liquid nitrogen, and stored at -70°C until isolation.
- Frozen bacterial pellets should not be thawed prior to beginning the protocol. Add the Lysozyme-containing TE Buffer directly to the frozen bacterial pellet (step **b** under lysate preparation).
- All centrifugation steps involving a microcentrifuge are carried out at $14,000 \times g$ ($\sim 12,000$ RPM) except where noted. All centrifugation steps are performed at room temperature.

Lysate Preparation

- a. Transfer 20 - 30 mL of urine to a 50 mL conical tube and centrifuge at $3,000 \times g$ for 5 minutes in a swinging bucket centrifuge to pellet the cells. Pour off the supernatant carefully so as not to disturb or dislodge the cell pellet.
- b. Resuspend the bacteria thoroughly in 100 μ L of the lysozyme-containing TE Buffer (prepared by user) by vortexing. Incubate at room temperature for 10 minutes.

Note: The length of the incubation step may be decreased to as little as 5 minutes if the urine bacteria of interest from which RNA is being isolated are known to be Gram-negative.

- c. Add 300 μ L of **Lysis Solution** and vortex vigorously for at least 10 seconds. Transfer the lysate to a microcentrifuge tube.
- d. Add 200 μ L of 95 – 100% ethanol (provided by the user) to the lysate. Mix by vortexing for 10 seconds.
- e. If any visible precipitates are present, spin the lysate for 1 minute in a benchtop microcentrifuge to pellet any debris. Otherwise, proceed directly to step **III A** without centrifugation.

III. Purifying Total RNA from Sample Lysate

A. Binding RNA to Column

- Assemble a spin column with one of the provided collection tubes.
- Apply the clarified lysate with the ethanol onto the column and centrifuge for 1 minute.
- Discard the flowthrough. Reassemble the spin column with its collection tube.

Optional Step: The SurePrep™ Urine Bacterial RNA Purification Kit isolates total RNA with minimal amounts of genomic DNA contamination. However, an optional **On-Column DNA Removal Protocol** is provided in Appendix A for maximum removal of residual DNA that may affect sensitive downstream applications. This optional step should be performed at this point in the protocol.

B. Column Wash

- Apply 500 μL of **Wash Solution** to the column and centrifuge for 1 minute.

Note: Ensure the entire wash solution has passed through into the collection tube by inspecting the column. If the entire wash volume has not passed, spin for an additional minute.

- Discard the flowthrough and reassemble the spin column with its collection tube.
- Wash column a second time by adding another 500 μL of **Wash Solution** and centrifuging for 2 minutes.
- Ensure that the column is dry. Spin for an additional minute, if necessary.
- Discard the collection tube with the flowthrough.

C. RNA Elution

- Place the column into a fresh 1.7 mL elution tube provided with the kit.
- Add 50 μL of **RNA Elution Buffer** to the column.
- Centrifuge for 2 minutes at **200 x g (~1,500 RPM)**, followed by a 1 minute spin at **14,000 x g (~12,000 RPM)**. Note the volume eluted from the column. If the entire 50 μL has not been eluted, spin the column at 14,000 x g (~12,000 RPM) for 1 additional minute.

Note: A smaller volume of RNA Elution Buffer may be used in order to obtain a more concentrated sample. A minimum volume of 20 μL is recommended

D. Assessing RNA Yield and Quality by UV Absorbance

The concentration and purity of an RNA solution can be determined by absorbance (A) measurements at 260 and 280 nm. A_{260} measurements are quantitative for relatively pure RNA preparations in microgram quantities. A_{260} readings cannot distinguish between DNA and RNA, however the ratio of A_{260}/A_{280} can be used as an indication of RNA purity. For example, contaminating proteins have a peak absorption at 280 nm that will reduce the A_{260}/A_{280} ratio.

- Determine RNA concentration by diluting an aliquot of the purified RNA solution in TE (10 mM Tris and 1 mM EDTA, pH 7.4). Measure absorbance of the diluted sample in a 1 mL cuvette using a traditional UV-VIS spectrophotometer

at 260 and 280 nm. The spectrophotometer should first be zeroed with the TE used to dilute the sample.

- b. An A_{260} of 1.0 is equivalent to 40 $\mu\text{g RNA/mL}$. Calculate the RNA concentration in $\mu\text{g/mL}$ as follows:

$$A_{260} \times \text{dilution factor} \times 40 = \mu\text{g RNA/mL}$$

- c. The ratio of the readings at 260 and 280 nm (A_{260}/A_{280}) provides an estimate of the RNA purity with respect to contaminants that absorb in the UV range such as protein. Ratios of 1.8 to 2.1 indicate highly purified preparations of RNA. Contaminants such as protein that absorb at 280 nm will lower this ratio. However, RNA solutions with a ratio lower than 1.8 may function well in downstream applications such as RT-PCR and Northern blotting.

E. Storage of RNA

The purified RNA sample may be stored at -20°C for a few days. It is recommended that samples be placed at -70°C for long term storage.

IV. Troubleshooting Guide

Problem	Possible Cause	Solution and Explanation
Poor RNA Recovery	Incomplete lysis of bacterial cells	Ensure that the appropriate amount of lysozyme-containing TE buffer and Lysis Solution are added to the bacterial pellet in order to completely lyse the cells.
	Column has become clogged	Do not exceed the recommended amount of 50 mL of urine. The amount of starting material may need to be decreased if the column shows clogging below the recommended levels. See also “Clogged Column” below.
	An alternative RNA elution solution was used	It is recommended that the RNA Elution Buffer supplied with this kit be used for maximum RNA recovery.
	Ethanol was not added to the lysate	Ensure that the appropriate amount of ethanol is added to the lysate before binding to the column.
	Ethanol was not added to the Wash Solution	Ensure that 50 mL of 95% ethanol is added to the supplied Wash Solution prior to use.
	There is very little or no bacteria in the urine	The expected amount of bacteria in a urine sample is very little. A healthy individual usually has < 10,000 CFU/mL, therefore it is possible that the urine sample has very little bacteria present. The isolated RNA may not be visible when resolved on a gel. In such cases, a larger input volume may be used. Alternatively, a more sensitive method such as BioAnalyzer or PCR amplification may be used for detection.
Clogged Column	Insufficient solubilization of cells	Ensure that the appropriate amount of lysozyme-containing TE buffer and Lysis Solution are added to the bacterial pellet in order to completely lyse the cells.
	Too many bacteria present in the urine	The urine sample that was applied to the column contained too many bacterial cells. Reduce the amount of urine used.
	Clarified lysate was not used for the binding step	If particulates are present in the lysate ensure that the sample is centrifuged and only the clarified lysate is used in subsequent steps.
	Centrifuge temperature too low	Ensure that the centrifuge remains at room temperature throughout the procedure. Temperatures below 20°C may cause precipitates to form that can cause the columns to clog.

Problem	Possible Cause	Solution and Explanation
RNA is Degraded	RNase contamination	RNases may be introduced during the use of the kit. Ensure proper procedures are followed when working with RNA. Please refer to “ <i>Working with RNA</i> ” at the beginning of this user guide.
	Procedure not performed quickly enough	In order to maintain the integrity of the RNA, it is important that the procedure be performed quickly.
	Improper storage of the purified RNA	For short term storage RNA samples may be stored at –20°C for a few days. It is recommended that samples be stored at –70°C for longer term storage.
	Lysozyme used may not be RNase-free	Ensure that the lysozyme being used with this kit is RNase-free in order to prevent possible problems with RNA degradation.
	The urine sample is old.	The use of fresh urine samples is recommended. Proteases and RNases may be present in the sample, and storing the sample for too long before RNA isolation increases the chances of recovering degraded RNA.
RNA does not perform well in downstream applications	RNA was not washed twice with the provided Wash Solution	Traces of salt from the binding step may remain in the sample if the column is not washed twice with Wash Solution . Salt may interfere with downstream applications, and thus must be washed from the column.
	Ethanol carryover	Ensure that the dry spin under the Column Wash procedure is performed, in order to remove traces of ethanol prior to elution. Ethanol is known to interfere with many downstream applications.
Genomic DNA contamination	Large amounts of starting material used	Perform RNase-free DNase I digestion on the RNA sample after elution to remove genomic DNA contamination (Appendix A).

V. Related Products

A. Additional RNA Purification Kits

Catalog #	Product Description
BP2800-50	SurePrep™ TrueTotal™ RNA Purification Kit
BP2801-25	SurePrep™ Small RNA Purification Kit
BP2802-50	SurePrep™ RNA/DNA/Protein Purification Kit
BP2803-50	SurePrep™ Urine Exfoliated Cell RNA Purification Kit
BP2804-50	SurePrep™ Urine Bacterial RNA Purification Kit
BP2805-50	SurePrep™ Nuclear Or Cytoplasmic RNA Purification Kit
BP2806-50	SurePrep™ RNA/Protein Purification Kit
BP2807-50	SurePrep™ Leukocyte RNA Purification Kit
BP2809-50	SurePrep™ RNA Cleanup and Concentration Kit

B. Other Fisher BioReagents Functionally Tested for RNA Research

BP2484-50	Water, Sterile (DEPC-treated) 50mL
BP2484-100	Water, Sterile (DEPC-treated) 100mL
BP561-1	Water, Sterile (RNA Grade) 1L
BP2483-100	EDTA 0.5 M (DEPC-treated) 100mL
BP2483-1	EDTA 0.5 M (DEPC-treated) 1L
BP2483-500	EDTA 0.5 M (DEPC-treated) 500mL
BP2810-50	RiboLadder™ 100b RNA Standard with loading buffers
BP2811-50	RiboLadder™ 1Kb RNA Standard with loading buffers
BP3224-5	Optizyme™ Ribonuclease Inhibitor (Human Placental) 10,000U
BP3224-1	Optizyme™ Ribonuclease Inhibitor (Human Placental) 2,500U
BP3225-5	Optizyme™ Ribonuclease Inhibitor (Porcine) 10,000U
BP3225-1	Optizyme™ Ribonuclease Inhibitor (Porcine) 2,500U
BP3222-5	Optizyme™ Ribonuclease Inhibitor (Recombinant) 10,000U
BP3222-1	Optizyme™ Ribonuclease Inhibitor (Recombinant) 2,500U
BP3226-1	Optizyme™ Recombinant DNase I (RNase-free) 1,000U
BP3226-2	Optizyme™ Recombinant DNase I (RNase-free) 2,000U
BP176-100	2-Mercaptoethanol 100g
BP535-1	Lysozyme, Egg White 1g
BP535-5	Lysozyme, Egg White 5g
BP535-10	Lysozyme, Egg White 10g
BP2476-100	Tris-EDTA, 1X Solution, pH 7.4 100ml
BP2476-500	Tris-EDTA, 1X Solution, pH 7.4 500ml
BP160-100	Agarose, Low EEO, Multipurpose 100g
BP1360-100	Agarose, Low Melting, <1kb RNA 100g
BP1356-100	Agarose, Broad Separation Range for RNA 100g
BP308-100	MOPS 100g
BP308-500	MOPS 500g

VI. Appendix A

Protocol for Optional On-Column DNA Removal

The SurePrep™ Urine Bacterial RNA Purification Kit isolates total RNA with minimal amounts of genomic DNA contamination. However, an optional protocol is provided below for maximum removal of residual DNA that may affect sensitive downstream applications. It is recommended that an RNase-free DNase I is used.

1. Prepare a working stock of 0.25 Kunitz unit/ μ L RNase-free DNase I solution according to the manufacturer's instructions. A 100 μ L aliquot is required for each column to be treated. Alternatively, dissolve or dilute stock DNase I in a reaction buffer (40 mM Tris pH 8.0, 10 mM MgCl₂ and 3 mM CaCl₂, made RNase-free) to give a final concentration of 0.25 Kunitz unit/ μ L.
2. Perform the procedure up to and including "**Binding RNA to Column**" (step III A).
3. Apply 400 μ L of **Wash Solution** to the column and centrifuge for 2 minutes. Discard the flowthrough. Reassemble the spin column with its collection tube.
4. Apply 100 μ L of the RNase-free DNase I solution prepared in step 1 above to the column. Centrifuge for 30 seconds at 200 x g (~1500 RPM). Alternatively, centrifuge for a 5 second pulse at 14, 000 x g (~12 000 RPM) if only a single speed centrifuge is available. Approximately one half of the DNase I solution will pass through the column.
5. Incubate the column assembly at 25-30°C for 15 minutes.
6. Without further centrifugation, proceed directly to "**Column Wash**" (step III B).

Technical Support

Telephone: 1-800-766-7000

E-mail: chem.techinfo@thermofisher.com

© 2007 Fisher BioReagents®

BP2007-04