



SurePrep™ Soil DNA Isolation Kit

Product Cat. # BP2815-50

Instruction Manual

| | | |
|-------------|--|----------|
| I. | Introduction | 2 |
| | A. Product Description | |
| | B. Overview of Procedure | |
| | C. Kit Specifications | |
| | D. Advantages | |
| | E. Kit Components | |
| | F. Storage Conditions and Product Stability | |
| | G. Precautions and Disclaimers | |
| | H. Customer-Supplied Reagents and Equipment | |
| II. | Set-Up and Preparation of Sample Lysate | 5 |
| | A. Equipment Preparation | |
| | B. Preparation of Lysate from Soil Samples | |
| III. | Purifying DNA from Soil Sample Lysates | 7 |
| | A. Binding DNA to Column | |
| | B. Column Wash | |
| | C. DNA Elution | |
| | D. Assessing DNA Yield by UV Absorbance | |
| | E. Storage of DNA Samples | |
| IV. | Troubleshooting Guide | 8 |
| V. | Related Products | 9 |
| | A. Additional RNA Purification Kits | |
| | B. Other Fisher BioReagents Functionally Tested for RNA Research | |

I. Introduction

A. Product Description

The SurePrep™ Soil DNA Isolation Kit provides a convenient and rapid method for the detection of microorganisms from soil samples. All types of soil samples can be processed with this kit, including common soil samples and difficult soil samples with high humic acid content such as compost and manure. The kit removes all traces of humic acid using the provided beads and a combination of chemical and physical homogenization and lysis. A simple and rapid spin column procedure is then used to further purify the DNA. Total genomic DNA can be isolated and purified from all the various microorganisms found in soil, such as bacteria, fungi and algae. The purified DNA is of the highest quality and is fully compatible with downstream PCR applications, since all humic acid substances and PCR inhibitors are removed during the isolation.

B. Overview of Procedure

Purification is based on spin column chromatography using Fisher's proprietary resin as the separation matrix. The process involves first measuring the soil sample into a provided 2 mL Screw Cap Tube, and adding an appropriate amount of the provided Beads. Next, Lysis Solution is added to the tube and vortexed briefly to mix. Lysis Additive is then added and the tube is vortexed for 5 minutes in order to efficiently and rapidly homogenize the sample, extract the DNA and remove all humic acids. The sample is then centrifuged, and the supernatant is transferred to a DNase-free microcentrifuge tube. Binding Solution is added, and the lysate is incubated for 5 minutes on ice. The lysate is then spun for 5 minutes to pellet any cell debris, the supernatant is collected, an equal volume of ethanol is added to the lysate and the solution is loaded onto a spin-column. Fisher's resin binds nucleic acids in a manner that depends on ionic concentrations, thus only the DNA will bind to the column while the proteins are removed in the flowthrough or retained on top of the resin. The bound DNA is then washed using the provided Wash Solutions, and the purified DNA is eluted using the Elution Buffer. The purified total DNA is free of all inhibitors, including humic acid, and can be used in sensitive downstream applications such as PCR.

C. Kit Specifications

| Kit Specifications | |
|---|---|
| Maximum Soil Input: Clay, loam and sand Feces and compost | 250 mg 100-150 mg |
| Type of Soil Processed | All types, including common soil, compost and manure |
| Column Binding Capacity | 50 µg |
| Maximum Column Loading Volume | 600 µL |
| Time to Complete 10 Purifications | 30 minutes |

D. Advantages

- Rapid and convenient method to detect microorganisms in soil samples
- Process all types of soil, including common soil, compost and manure
- Remove all humic acid from DNA samples
- Fast and easy processing using a rapid spin-column format
- Isolate high quality total DNA from a variety of microorganisms including bacteria, fungi and algae

E. Kit Components

| Component | Catalog # BP2815-50 (50 preps) |
|------------------------|--------------------------------|
| Lysis Solution | 45 mL |
| Lysis Additive | 6 mL |
| Binding Solution | 6 mL |
| Wash Solution I | 30 mL |
| Wash Solution II | 18 mL |
| Elution Buffer | 6 mL |
| Beads | 30 g |
| Screw Cap Tubes (2 mL) | 50 |
| Spin Columns | 50 |
| Collection Tubes | 50 |
| Elution tubes (1.7 mL) | 50 |
| Product Insert | 1 |

F. Storage Conditions and Product Stability

All kit solutions should be kept tightly sealed and stored at room temperature. These reagents should remain stable for at least 2 years in their unopened containers.

G. Precautions and Disclaimers

This kit is designed for research purposes only. It is not intended for human or diagnostic use.

Ensure that a suitable lab coat, disposable gloves and protective goggles are worn when working with chemicals. For more information, please consult the appropriate Material Safety Data Sheets (MSDS). The MSDS can be requested through our Customer Service Department.

H. Customer-Supplied Reagents and Equipment


You must have the following in order to use the SurePrep™ Soil DNA Isolation Kit:

- Benchtop microcentrifuge
- DNase-free microcentrifuge tubes
- Flat bed vortex or bead beater equipment
- 95-100% ethanol
- 70% ethanol

Flowchart

*Procedure for Purifying Total DNA Using
the SurePrep Soil DNA Isolation Kit*

Add beads, soil sample, Lysis Solution and Lysis Additive to screw cap tube



II. Set-Up and Preparation of Sample Lysate

A. Equipment Preparation

All centrifugation steps are carried out in a benchtop microcentrifuge. Various speeds are required for different steps, so please check your microcentrifuge specifications to ensure that it is capable of the proper speeds. All centrifugation steps are performed at room temperature. The correct rpm can be calculated using the formula:

$$RPM = \sqrt{\frac{RCF}{(1.118 \times 10^{-5}) (r)}}$$

where *RCF* = required gravitational acceleration (relative centrifugal force in units of g); *r* = radius of the rotor in cm; and *RPM* = the number of revolutions per minute required to achieve the necessary *g*-force.

If you do not own a variable speed microcentrifuge consider purchasing Thermo Scientific's Sorvall Legend benchtop model that offers power, safety and convenience.

- Choice of 17,000 or 21,000 x g (for RNA purification using SurePrep kits, the microcentrifuge with 17,000 x g is sufficient)
- Holds 36 x 0.5 mL microtubes, 24 x 2 mL tubes or 8 x 8 PCR
- Unique ClickSeal™ bio-containment rotor lid for safe processing of infectious specimens
- Fast acceleration and deceleration speeds up your protocols
- Broad range of rotors supports virtually any application
- Intuitive controls and vivid display
- Highly resistant materials allow vigorous cleaning and autoclaving

Sorvall Legend Micro Centrifuges

Technical Specifications

| | Sorvall Legend Micro 17 & 17R | Sorvall Legend Micro 21 & 21R |
|-----------------|--|---|
| Max g-force: | 17,000 | 21,100 |
| Max RPM: | 13,300 | 14,800 |
| Noise level: | <55 dBA | <56 dBA |
| Time set range: | 1 min - 99 min; 1 min increments | 1 min - 99 min; 1 min increments |
| Temp set range: | Set from -9 °C to +40 °C; per 1 °C increment | Set from -9 °C to +40 °C per 1 °C increment |

Ordering Information

| | Cat. No. | Cat. No. |
|---|--------------|------------|
| Sorvall Legend Micro 17/17R | 230V 50/60Hz | 120V 60 Hz |
| Sorvall Legend Micro 17, includes 24 x 1.5/2.0 mL rotor with ClickSeal bio-containment lid | 75002430 | 75002431 |
| Sorvall Legend Micro 17R, includes 24 x 1.5/2.0 mL rotor with ClickSeal bio-containment lid | 75002440 | 75002441 |
| Sorvall Legend Micro 21/21R | 230V 50/60Hz | 120V 60 Hz |
| Sorvall Legend Micro 21, incl. 24 x 1.5/2.0 mL rotor with ClickSeal bio-containment lid | 75002435 | 75002436 |
| Sorvall LegendMicro 21R, incl. 24 x 1.5/2.0 mL rotor with ClickSeal bio-containment lid | 75002445 | 75002446 |

For detailed product specifications, information on additional rotors, lids and adapters visit www.thermo.com

B. Preparation of Lysate from Soil Samples

Notes Prior to Use

- All centrifugation steps are carried out in a benchtop microcentrifuge at **14,000 x g (~12,000 RPM)** except where noted. All centrifugation steps are performed at room temperature.
- A variable speed centrifuge should be used for maximum kit performance. If a variable speed centrifuge is not available a fixed speed centrifuge can be used, however reduced yields may be observed.
- Ensure that all solutions are at room temperature prior to use.
- Prepare a working concentration of **Wash Solution II** by adding 42 mL of 95 - 100 % ethanol (provided by the user) to the supplied bottle containing the concentrated **Wash Solution II**. This will give a final volume of 60 mL. The label on the bottle has a box that may be checked to indicate that the ethanol has been added.
- Recommended soil input varies depending on the soil type. For clay, loam and sand the recommended soil input is 250 mg. For fecal samples and compost, it is recommended that 100-150 mg of sample be used.

Lysate Preparation

- a. Weigh up to 250 mg of soil sample (please see Notes Prior to Use) into a 2 mL **Screw Cap Tube** provided with the kit.

Note: In case of a wet soil sample, transfer the sample to a clean 1.7 mL microcentrifuge tube and centrifuge for 30 seconds at **14,000 x g (~12,000 RPM)**. Remove the water carefully using a pipette, and weigh the remaining soil.

- b. Add a 2:1 ratio of the provided **Beads** into the screw cap tube with the soil (500 mg of glass beads is added to 250 mg of soil).
- c. Add 700 μ L of **Lysis Solution** to the tube. Vortex briefly to mix soil and Lysis Solution.
- d. Add 100 μ L of **Lysis Additive** and vortex briefly.
- e. Secure tube horizontally on a flat-bed vortex pad with tape, or secure the tube in any commercially available bead beater equipment (e.g. Scientific Industries' Disruptor Genie™). Vortex for 5 minutes at maximum speed.
- f. Centrifuge the tube for 1 minute at **14,000 x g (~12,000 RPM)**.
- g. Transfer up to 450 μ L of supernatant to a DNase-free microcentrifuge tube (not provided).
- h. Add 100 μ L of **Binding Solution**, mix by inverting the tube a few times, and incubate for 5 minutes on ice.
- i. Spin the lysate for 1 minute to pellet any protein and soil particles.
- j. Using a pipette, transfer up to 450 μ L of supernatant into a DNase-free microcentrifuge tube (not provided).

Note: Avoid any contact with the pellet when collecting the supernatant. Also, depending on the soil type, some residue may be present on top of the supernatant. It is important to avoid collection of this residue while collecting the supernatant.

- k. Add an equal volume of 70% ethanol (provided by the user) to the lysate collected above (100 μ L of ethanol is added to every 100 μ L of lysate). Vortex to mix. Proceed to **step III. A Binding DNA to Column.**

III. Purifying DNA from Soil Sample Lysates

A. Binding DNA to Column

- a. Assemble a **spin column** with one of the provided **collection tubes**.
- b. Apply up to 600 μL of the clarified lysate with ethanol onto the column and centrifuge for 1 minute at **14,000 \times g (~12,000 RPM)**. Discard the flowthrough and reassemble the spin column with the collection tube.

Note: Ensure the entire lysate volume has passed through into the collection tube by inspecting the column. If the entire lysate volume has not passed, spin for an additional minute.

- c. Depending on your lysate volume, repeat step **b** if necessary.

B. Column Wash

- a. Apply 500 μL of **Wash Solution I** to the column and centrifuge for 1 minute.

Note: Ensure the entire wash solution has passed through into the collection tube by inspecting the column. If the entire wash volume has not passed, spin for an additional minute.

- b. Discard the flowthrough and reassemble the spin column with its collection tube.
- c. Apply 500 μL of **Wash Solution II** to the column and centrifuge for 1 minute.
- d. Discard the flowthrough and reassemble the spin column with its collection tube.
- e. Repeat **c** and **d**.
- f. Spin the column for 2 minutes in order to thoroughly dry the resin. Discard the collection tube.

C. DNA Elution

- a. Place the column into a fresh 1.7 mL **Elution Tube** provided with the kit.
- b. Add 50 μL of **Elution Buffer** to the column.
- c. Centrifuge for 2 minutes at **200 \times g (~1,500 RPM)**, followed by a one minute spin at **14,000 \times g (~12,000 RPM)**. Note the volume eluted from the column. If the entire volume has not been eluted, spin the column at 14,000 \times g (~12,000 RPM) for one additional minute.
- d. **(Optional)** An additional elution may be performed if desired by repeating steps **b** and **c** using 50 μL of **Elution Buffer**. The total yield can be improved by an additional 20-30% when this second elution is performed.

D. Assessing DNA Yield by UV Absorbance

Spectrophotometric measurement of the amount of UV irradiation absorbed by DNA is simple and accurate when the DNA sample is pure.

- a. Determine DNA concentration by making an appropriate dilution of the purified DNA solution in TE (10 mM Tris and 1 mM EDTA, pH 7.4). Measure absorbance of the diluted sample in a 1 mL cuvette using a traditional UV-VIS spectrophotometer at 260 and 280 nm. The spectrophotometer should first be zeroed with the TE used to dilute the sample.
- b. An A_{260} of 1.0 is equivalent to 50 μg double-stranded DNA per mL. Calculate the DNA concentration in $\mu\text{g}/\text{mL}$ as follows:

$$A_{260} \times \text{dilution factor} \times 50 = \mu\text{g DNA/mL}$$

- c. The ratio of the readings at 260 and 280 nm (A_{260}/A_{280}) provides an estimate of the DNA purity with respect to contaminants that absorb in the UV range such as protein. Ratios of 1.8 to 2.0 indicate highly purified preparations of DNA. Contaminants such as protein that absorb at 280 nm will lower this ratio.

E. Storage of DNA Samples

The purified genomic DNA can be stored at 2-8°C for a few days. For longer term storage, -20°C is recommended.

IV. Troubleshooting Guide

| Problem | Possible Cause | Solution and Explanation |
|--|---|--|
| Poor DNA Recovery | Homogenization was incomplete | Depending on the type of soil, further vortexing with the flat bed vortex or bead beater equipment may be required. However, it is not recommended to increase the vortex time to longer than 10 minutes at maximum speed. Also, ensure that the maximum recommended soil input is not exceeded, as this may also cause incomplete homogenization. |
| | An alternative elution buffer was used | It is recommended that the Elution Buffer supplied with this kit be used for maximum DNA recovery. |
| | Lysis Additive was not added to the lysate | Ensure that the provided Lysis Additive is added to separate humic acid and increase DNA yield. |
| | Binding Solution was not added to the lysate | Ensure that the Binding Solution is added to the lysate and that it is incubated on ice for 5 minutes prior to spinning down the lysate. |
| | Ethanol was not added to the lysate | Ensure that an equal amount of ethanol is added to the lysate before binding to the column. |
| | Ethanol was not added to the Wash Solution II | Ensure that 42 mL of 95 - 100% ethanol is added to the supplied Wash Solution II prior to use. |
| DNA does not perform well in downstream applications | Eluted DNA sample is brown | Ensure that the Lysis Additive is added. Also, avoid any contact with the pellet or surface residue when collecting the supernatant during sample preparation. |
| | Too much soil input was processed | Ensure that the maximum input of soil is not exceeded, as this may also cause incomplete homogenization. |
| | Lysis Additive was not added to the lysate | Ensure that the provided Lysis Additive is added to the lysate. Also, an incubation can be performed at 65°C for 10 minutes after addition of the Lysis Additive and prior to vortexing to maximize DNA recovery. |

| Problem | Possible Cause | Solution and Explanation |
|--|---|--|
| DNA does not perform well in downstream applications | DNA was not washed three times with the provided Wash Solutions | Traces of salt from the binding step may remain in the sample if the column is not washed three times with the provided Wash Solutions . Salt may interfere with downstream applications, and therefore must be washed from the column. |
| | Ethanol carryover | Ensure that the dry spin under the Column Wash procedure is performed, in order to remove traces of ethanol prior to elution. Ethanol is known to interfere with many downstream applications. |
| | PCR reaction conditions need to be optimized | Take steps to optimize the PCR conditions being used, including varying the amount of template, changing the source of <i>Taq</i> polymerase, verifying the primer design and adjusting the primer annealing conditions. |

V. Related Products

A. Additional RNA Purification Kits

| Catalog # | Product Description |
|------------------|---|
| BP2800-50 | SurePrep™ TrueTotal™ RNA Purification Kit |
| BP2801-25 | SurePrep™ Small RNA Purification Kit |
| BP2802-50 | SurePrep™ RNA/DNA/Protein Purification Kit |
| BP2803-50 | SurePrep™ Urine Exfoliated Cell RNA Purification Kit |
| BP2804-50 | SurePrep™ Urine Bacterial RNA Purification Kit |
| BP2805-50 | SurePrep™ Nuclear Or Cytoplasmic RNA Purification Kit |
| BP2806-50 | SurePrep™ RNA/Protein Purification Kit |
| BP2807-50 | SurePrep™ Leukocyte RNA Purification Kit |
| BP2809-50 | SurePrep™ RNA Cleanup and Concentration Kit |
| BP2814-25 | SurePrep™ Water RNA/DNA Purification Kit |
| BP2815-50 | SurePrep™ Soil DNA Isolation Kit |
| BP2816-50 | SurePrep™ FFPE RNA Isolation Kit |
| BP2817-50 | SurePrep™ Plant/Fungi Total RNA Purification Kit |

B. Other Fisher BioReagents Functionally Tested for RNA Research

| | |
|------------|--|
| BP2484-50 | Water, Sterile (DEPC-treated) 50mL |
| BP2484-100 | Water, Sterile (DEPC-treated) 100mL |
| BP561-1 | Water, Sterile (RNA Grade) 1L |
| BP2483-100 | EDTA 0.5 M (DEPC-treated) 100mL |
| BP2483-1 | EDTA 0.5 M (DEPC-treated) 1L |
| BP2483-500 | EDTA 0.5 M (DEPC-treated) 500mL |
| BP2810-50 | RiboLadder™ 100b RNA Standard with loading buffers |
| BP2811-50 | RiboLadder™ 1Kb RNA Standard with loading buffers |

| | |
|------------|--|
| BP3224-5 | Optizyme™ Ribonuclease Inhibitor (Human Placental) 10,000U |
| BP3224-1 | Optizyme™ Ribonuclease Inhibitor (Human Placental) 2,500U |
| BP3225-5 | Optizyme™ Ribonuclease Inhibitor (Porcine) 10,000U |
| BP3225-1 | Optizyme™ Ribonuclease Inhibitor (Porcine) 2,500U |
| BP3222-5 | Optizyme™ Ribonuclease Inhibitor (Recombinant) 10,000U |
| BP3222-1 | Optizyme™ Ribonuclease Inhibitor (Recombinant) 2,500U |
| BP3226-1 | Optizyme™ Recombinant DNase I (RNase-free) 1,000U |
| BP3226-2 | Optizyme™ Recombinant DNase I (RNase-free) 2,000U |
| BP176-100 | 2-Mercaptoethanol 100g |
| BP535-1 | Lysozyme, Egg White 1g |
| BP535-5 | Lysozyme, Egg White 5g |
| BP535-10 | Lysozyme, Egg White 10g |
| BP2476-100 | Tris-EDTA, 1X Solution, pH 7.4 100ml |
| BP2476-500 | Tris-EDTA, 1X Solution, pH 7.4 500ml |
| BP160-100 | Agarose, Low EEO, Multipurpose 100g |
| BP1360-100 | Agarose, Low Melting, <1kb RNA 100g |
| BP1356-100 | Agarose, Broad Separation Range for RNA 100g |
| BP1302-10 | Ethidium Bromide, 1% Solution 10ml |
| BP2900-500 | MOPS 10X Solution 500ml |
| BP2900-1 | MOPS 10X Solution 1L |

Fisher BioReagents Technical Support

Telephone: 1-800-766-7000

E-mail: chem.techinfo@thermofisher.com

© 2009 Fisher BioReagents®

Rev: New, 3/09