



Inorganic Salts Guide

Aluminum | Ammonium | Barium | Calcium | Lithium | Magnesium | Potassium | Sodium

Reliability.
Purity.
Certainty.



Table of Contents

Aluminum.....	3	Cesium.....	6	Rubidium.....	9
Ammonium.....	3-4	Lithium.....	6	Sodium.....	9-15
Barium.....	4-5	Magnesium.....	6	Strontium.....	15
Calcium.....	5-6	Potassium.....	7-9		

The Fisher Chemical Brand...



has a wide selection of inorganic salts that can be utilized in a number of diverse applications. These are offered in various purity grades to meet both your research and manufacturing process needs. Inorganic salts come in a variety of forms, including powder, crystals, granules, pellets and lumps. Grades of these inorganic salts include select products meeting ACS, USP/NF and FCC requirements along with Certified, Reagent and other specialty grades. All Fisher Chemical products are manufactured in an ISO 9001:2008 certified environment and undergo rigorous quality assurance and testing procedures ensuring excellent lot-to-lot consistency.

Table 1: Purity Grade Classification for Inorganic Salts

Grade	Definition
Multicompendia (USP/NF/FCC/EP/BP/JP)	Reagent chemicals that meet or surpass specifications of the United States Pharmacopeia (USP), the National Formulary (NF), the Food Chemicals Codex (FCC), the European Pharmacopeia (EP), the British Pharmacopeia (BP) and/or the Japanese Pharmacopeia (JP).
Certified ACS	Reagent chemicals that meet or exceed the latest ACS specifications. Actual lot analysis on label.
Certified	Reagent chemicals for which the purity standard is established by Fisher Chemical. Purity is guaranteed to meet the published maximum limits of impurities.
HPLC	Compounds manufactured specifically for use with HPLC instruments. Meets all ACS specifications. Submicron filtered.
Laboratory and Technical	Chemicals of reasonable quality and purity for use in procedures where no official standards are required.

Specialized Chemical Services

Our Specialized Chemical Services provide tailored chemical solutions and services required in research, scale up/development and production. We offer bulk and semi bulk quantities, solvent blending, additional testing, packaging and/or special instructions. Our robust Global Quality System provides cGMP products and services including supply chain integrity and traceability.



Ammonium

Product Description	Grade	Form	Size	Cat. No.
Aluminum				
Aluminum Ammonium Sulfate Dodecahydrate	Certified	Crystals	500g	A567500
Aluminum Chloride Hexahydrate	USP	Crystals	500g	A576500
Aluminum Chloride Hexahydrate	USP	Crystals	2.5kg	A576212
Aluminum Chloride Hexahydrate	Certified	Crystals	500g	A573500
Aluminum Chloride Hexahydrate	Certified	Crystals	2.5kg	A573212
Aluminum Hydroxide Hydrate	Certified	Powder	500g	A581500
Aluminum Nitrate Nonahydrate	ACS Certified	Crystals	250g	A586250
Aluminum Nitrate Nonahydrate	ACS Certified	Crystals	500g	A586500
Aluminum Nitrate Nonahydrate	ACS Certified	Crystals	3kg	A5863
Aluminum Nitrate Nonahydrate	ACS Certified	Crystals	10kg	A58610
Aluminum Potassium Sulfate Dodecahydrate	USP/FCC	Powder	500g	A605500
Aluminum Potassium Sulfate Dodecahydrate	USP/FCC	Powder	2.5kg	A605212
Aluminum Potassium Sulfate Dodecahydrate	ACS Certified	Crystals	500g	A601500
Aluminum Potassium Sulfate Dodecahydrate	ACS Certified	Crystals	3kg	A6013
Aluminum Sulfate Hydrate	Technical		3kg	A6113
Aluminum Sulfate Hydrate	ACS Certified	Crystals	500g	A613500
Ammonium				
Ammonium Acetate	ACS Certified	Crystals	500g	A637500
Ammonium Acetate	ACS Certified	Crystals	3kg	A6373
Ammonium Acetate	ACS Certified	Crystals	10kg	A63710
Ammonium Acetate	HPLC	Crystals	500g	A639500
Ammonium Bicarbonate	Certified	Powder	500g	A643500
Ammonium Bicarbonate	Certified	Powder	50kg	A64350
Ammonium Bifluoride	Technical	Crystals	500g	A664500
Ammonium Bifluoride	Technical	Crystals	3kg	A6643
Ammonium Bromide	ACS Certified	Powder	500g	A650500
Ammonium Carbonate	NF	Lumps	500g	A657500
Ammonium Carbonate	NF	Lumps	2.5kg	A657212
Ammonium Carbonate	NF	Lumps	12kg	A65712
Ammonium Carbonate	ACS Certified	Lumps	500g	A656500
Ammonium Carbonate	ACS Certified	Lumps	3kg	A6563
Ammonium Carbonate	HPLC	Lumps	500g	A651500
Ammonium Carbonate	Laboratory		3kg	A6523
Ammonium Chloride	USP/FCC/EP		500g	A687500
Ammonium Chloride	USP/FCC/EP		2.5kg	A687212
Ammonium Chloride	USP/FCC/EP		10kg	A68710
Ammonium Chloride	USP/FCC	Granular	500g	A649500
Ammonium Chloride	USP/FCC	Granular	3kg	A6493
Ammonium Chloride	ACS Certified	Crystals	500g	A661500
Ammonium Chloride	ACS Certified	Crystals	3kg	A6613
Ammonium Chloride	ACS Certified	Crystals	10kg	A66110
Ammonium Citrate Dibasic	ACS Certified	Crystals	500g	A663500
Ammonium Citrate Dibasic	ACS Certified	Crystals	2.5kg	A663212
Ammonium Dichromate	Certified	Crystals	100g	A644100

Ammonium

Product Description	Grade	Form	Size	Cat. No.
Ammonium (Continued)				
Ammonium Dichromate	Certified	Crystals	500g	A644500
Ammonium Dichromate	Certified	Crystals	2.5kg	A644212
Ammonium Ferric Citrate	Laboratory		500g	I72500
Ammonium Fluoride	ACS Certified	Crystals	100g	A665100
Ammonium Fluoride	ACS Certified	Crystals	500g	A665500
Ammonium Formate	Certified	Crystals	500g	A666500
Ammonium Hydroxide	ACS Certified		385 lb.	A669385LB
Ammonium Hydroxide	ACS Certified		2.5L	A669212
Ammonium Hydroxide, Safecote	ACS Certified		6 x 2.5L	A669S212
Ammonium Hydroxide, Safecote	ACS Certified		500mL	A669S500
Ammonium Hydroxide	ACS Certified		500mL	A669500
Ammonium Hydroxide	ACS Certified		6 x 2.5L	A669C212
Ammonium Hydroxide	ACS Certified		6 x 5L	A669612 GAL
Ammonium Metavanadate	ACS Certified		500g	A713500
Ammonium Molybdate Tetrahydrate	ACS Certified	Crystals	500g	A674500
Ammonium Molybdate Tetrahydrate	ACS Certified	Crystals	3kg	A6743
Ammonium Molybdate Tetrahydrate	ACS Certified	Crystals	10kg	A67410
Ammonium Nitrate	ACS Certified	Granular	500g	A676500
Ammonium Nitrate	ACS Certified	Granular	2.5kg	A676212
Ammonium Oxalate Monohydrate	ACS Certified	Crystals	500g	A679500
Ammonium Peroxydisulfate	ACS Certified	Crystals	500g	A682500
Ammonium Peroxydisulfate	ACS Certified	Crystals	3kg	A6823
Ammonium Peroxydisulfate	ACS Certified	Crystals	10kg	A68210
Ammonium Phosphate Dibasic	ACS Certified	Crystals	500g	A686500
Ammonium Phosphate Monobasic	ACS Certified	Crystals	500g	A684500
Ammonium Phosphate Monobasic	ACS Certified	Crystals	3kg	A6843
Ammonium Phosphate Monobasic	HPLC	Crystals	500g	A685500
Ammonium Pyrrolidine Dithiocarbamate	Laboratory		10g	A18210
Ammonium Sulfate	ACS Certified	Granular	500g	A702500
Ammonium Sulfate	ACS Certified	Granular	3kg	A7023
Ammonium Sulfate	ACS Certified	Granular	10kg	A70210
Ammonium Sulfate	Laboratory	Crystals	500g	A938500
Ammonium Sulfate	Laboratory		3kg	A7013
Ammonium Sulfate	Laboratory		50kg	A70150
Ammonium Thiocyanate	ACS Certified	Crystals	500g	A709500
Ammonium Thiocyanate	ACS Certified	Crystals	3kg	A7093
Ammonium Vanadate	Laboratory		500g	A714500
Barium				
Barium Acetate	ACS Certified	Crystals	250g	B24250
Barium Acetate	ACS Certified	Crystals	500g	B24500
Barium Carbonate	ACS Certified	Powder	100g	B30100
Barium Carbonate	ACS Certified	Powder	500g	B30500
Barium Carbonate	ACS Certified	Powder	3kg	B303
Barium Chloride Anhydrous	Laboratory		100g	B31100C

Calcium

Product Description	Grade	Form	Size	Cat. No.
Barium Chloride Anhydrous	Laboratory		500g	B31500
Barium Chloride Dihydrate	Technical	Crystals	10kg	B3310
Barium Chloride Dihydrate	Technical	Crystals	50kg	B3350
Barium Chloride Dihydrate	ACS Certified	Crystals	100g	B34100
Barium Chloride Dihydrate	ACS Certified	Crystals	500g	B34500
Barium Chloride Dihydrate	ACS Certified	Crystals	3kg	B343
Barium Chloride Dihydrate	ACS Certified	Crystals	10kg	B3410
Barium Diphenylamine Sulfonate	Certified	Powder	25g	B3825
Barium Hydroxide Anhydrous	Certified	Powder	250g	B47250
Barium Hydroxide Anhydrous	Certified	Powder	500g	B47500
Barium Hydroxide Octahydrate	ACS Certified	Crystals	500g	B46500
Barium Nitrate	ACS Certified	Crystals	500g	B53500
Barium Sulfate	USP	Powder	500g	B75500
Barium Sulfate	USP	Powder	2.5kg	B75212
Barium Sulfate	USP	Powder	12kg	B7512
Barium Sulfate	Certified	Powder	500g	B68500
Calcium				
Calcium Acetate Monohydrate	Certified	Powder	500g	C46500
Calcium Carbonate	USP/FCC	Powder	3kg	C633
Calcium Carbonate	USP/FCC	Powder	10kg	C6310
Calcium Carbonate	ACS Certified	Powder	500g	C64500
Calcium Carbonate	ACS Certified	Powder	3kg	C643
Calcium Chloride	Certified	Desiccant, 20 Mesh and Finer	500g	C77500
Calcium Chloride	Certified	Desiccant, 20 Mesh and Finer	2.5kg	C77212
Calcium Chloride Anhydrous	Technical	Pellets, 4-20 Mesh, for Desiccators	500g	C614500
Calcium Chloride Anhydrous	Technical	Pellets, 4-20 Mesh, for Desiccators	3kg	C6143
Calcium Chloride Anhydrous	Technical	Pellets, 4-20 Mesh, for Desiccators	10kg	C61410
Calcium Chloride Dihydrate	USP/FCC/EP		500g	C69500
Calcium Chloride Dihydrate	USP/FCC	Fine, White Granules	500g	C70500
Calcium Chloride Dihydrate	ACS Certified		500g	C79500
Calcium Chloride Dihydrate	ACS Certified		3kg	C793
Calcium Chloride Hexahydrate	Certified	Crystals	500g	C78500
Calcium Fluoride	Certified	Powder	250g	C89250
Calcium Fluoride	Certified	Powder	500g	C89500
Calcium Hydride	Laboratory		100g	C94100
Calcium Hydroxide	USP/FCC	Powder	500g	C88500
Calcium Hydroxide	Certified	Powder	500g	C97500
Calcium Hydroxide	Certified	Powder	3kg	C973
Calcium Hydroxide	Certified	Powder	10kg	C9710
Calcium Hypochlorite	Certified	Granular	250g	C100250
Calcium Hypochlorite	Certified	Granular	500g	C100500
Calcium Lactate Pentahydrate	USP		500g	C110500
Calcium Nitrate Tetrahydrate	ACS Certified	Crystals	500g	C109500
Calcium Nitrate Tetrahydrate	ACS Certified	Crystals	3kg	C1093
Calcium Phosphate Dibasic	Laboratory	Reagent	500g	C135500

Calcium

Product Description	Grade	Form	Size	Cat. No.
Calcium (Continued)				
Calcium Phosphate Dibasic Dihydrate	USP/FCC	Powder	500g	C129500
Calcium Phosphate Dibasic Dihydrate	USP/FCC	Powder	12kg	C12912
Calcium Phosphate Dibasic Dihydrate	USP/FCC/EP/BP	Powder	12kg	C13012
Calcium Phosphate Tribasic	FCC/NF	Powder	500g	C133500
Calcium Phosphate Tribasic	FCC/NF	Powder	2.5kg	C133212
Calcium Stearate	Technical	Powder	3kg	C1363
Calcium Sulfate Hemihydrate	Laboratory	Calcined Powder	10kg	C13810
Cesium				
Cesium Chloride	Certified	Crystals	25g	C2425
Lithium				
Lithium Bromide	Laboratory	Fine Powder	500g	L117500
Lithium Carbonate	ACS Certified	Powder	500g	L119500
Lithium Chloride	Certified		100g	L121100
Lithium Chloride	Certified		500g	L121500
Lithium Chloride	Laboratory		500g	L120500
Lithium Fluoride	Certified	Powder	500g	L124500
Lithium Hydroxide Anhydrous	Laboratory	Powder	500g	L128500
Lithium Hydroxide Monohydrate	Laboratory	Crystals	500g	L127500
Lithium Nitrate	Certified	Powder	500g	L132500
Magnesium				
Magnesium Acetate Tetrahydrate	ACS Certified	Crystals	500g	M13500
Magnesium Carbonate	USP/FCC	Powder	500g	M27500
Magnesium Carbonate	USP/FCC	Powder	2.5kg	M27212
Magnesium Carbonate Anhydrous	Laboratory	Powder	3kg	M263
Magnesium Carbonate Basic Hydrate	Certified	Powder	100g	M29100
Magnesium Carbonate Basic Hydrate	Certified	Powder	500g	M29500
Magnesium Chloride Hexahydrate	USP/FCC	Crystals	500g	M35500
Magnesium Chloride Hexahydrate	USP/FCC	Crystals	2.5kg	M35212
Magnesium Chloride Hexahydrate	USP/FCC/EP/BP	Crystals	2.5kg	M36212
Magnesium Chloride Hexahydrate	USP/FCC	Crystals	12kg	M3512
Magnesium Chloride Hexahydrate	ACS Certified	Crystals	500g	M33500
Magnesium Chloride Hexahydrate	ACS Certified	Crystals	3kg	M333
Magnesium Hydroxide	USP/FCC	Powder	500g	M342500
Magnesium Nitrate Hexahydrate	ACS Certified	Crystals	500g	M46500
Magnesium Nitrate Hexahydrate	ACS Certified	Crystals	2.5kg	M46212
Magnesium Perchlorate Anhydrous	Certified		2.5kg	M54212
Magnesium Stearate	NF	Impalpable Powder	500g	M58500
Magnesium Stearate	NF	Impalpable Powder	2.5kg	M58212
Magnesium Stearate	NF	Impalpable Powder	50kg	M5850
Magnesium Sulfate Anhydrous	Certified	Powder	500g	M65500
Magnesium Sulfate Anhydrous	Certified	Powder	3kg	M653
Magnesium Sulfate Anhydrous	Certified	Powder	100 lb.	M65100LB
Magnesium Sulfate Heptahydrate	USP/FCC/EP		500g	M80500
Magnesium Sulfate Heptahydrate	USP/FCC/EP		2.5kg	M80212

Potassium

Product Description	Grade	Form	Size	Cat. No.
Magnesium Sulfate Heptahydrate	USP/FCC		3kg	M673
Magnesium Sulfate Heptahydrate	USP/FCC		50kg	M6750
Magnesium Sulfate Heptahydrate	ACS Certified	Crystals	500g	M63500
Magnesium Sulfate Heptahydrate	ACS Certified	Crystals	3kg	M633
Magnesium Sulfate Heptahydrate	ACS Certified	Crystals	50kg	M6350
Potassium				
Potassium Acetate	USP/EP/BP	Crystals	2.5kg	P181212
Potassium Acetate	USP	Crystals	3kg	P1783
Potassium Acetate	USP/EP/BP	Crystals	12kg	P18112
Potassium Acetate	ACS Certified	Crystals	500g	P171500
Potassium Acetate	ACS Certified	Crystals	3kg	P17113
Potassium Bicarbonate	USP/FCC	Powder	500g	P235500
Potassium Bicarbonate	USP/FCC	Powder	2.5kg	P235212
Potassium Bicarbonate	USP/FCC	Powder	12kg	P23512
Potassium Bicarbonate	ACS Certified	Crystals	500g	P184500
Potassium Biiodate	Laboratory		100g	P190100
Potassium Bisulfate	Certified	Crystals	500g	P193500
Potassium Bromate	Certified	With Magnesium Carbonate; Granular	250g	P207250
Potassium Bromide	Technical	Infrared	25g	P22725
Potassium Bromide	ACS Certified	Beads or Powder	500g	P205500
Potassium Bromide	ACS Certified	Beads or Powder	3kg	P2053
Potassium Carbonate Anhydrous	ACS Certified	Granular Powder	500g	P208500
Potassium Carbonate Anhydrous	ACS Certified	Granular Powder	3kg	P2083
Potassium Carbonate Anhydrous	ACS Certified	Granular Powder	250 lb.	P208250LB
Potassium Carbonate Sesquihydrate	ACS Certified	Powder	500g	P179500
Potassium Chlorate	ACS Certified	Crystals	100g	P212100
Potassium Chlorate	ACS Certified	Crystals	500g	P212500
Potassium Chlorate	Laboratory	Granular	500g	P210500
Potassium Chloride	EP/BP/USP/FCC	Crystals	500g	P333500
Potassium Chloride	USP/FCC	Crystals	500g	P330500
Potassium Chloride	USP/FCC/EP/BP/JP	Crystals	2.5kg	P335212
Potassium Chloride	EP/BP/USP/FCC	Crystals	3kg	P3333
Potassium Chloride	USP/FCC	Crystals	3kg	P3303
Potassium Chloride	USP/FCC/EP/BP/JP	Crystals	12kg	P33512
Potassium Chloride	EP/BP/USP/FCC	Crystals	250 lb.	P333250LB
Potassium Chloride	USP/FCC	Crystals	250 lb.	P330250LB
Potassium Chloride	ACS Certified	Crystals	500g	P217500
Potassium Chloride	ACS Certified	Crystals	3kg	P2173
Potassium Chloride	ACS Certified	Crystals	10kg	P21710
Potassium Chloride	ACS Certified	Crystals	250 lb.	P217250LB
Potassium Chromate	ACS Certified	Granular	100g	P220100
Potassium Chromate	ACS Certified	Granular	500g	P220500
Potassium Chromate	ACS Certified	Granular	3kg	P2203
Potassium Citrate Monohydrate	USP/FCC	Granular	500g	P218500
Potassium Citrate Monohydrate	USP/FCC	Granular	2.5kg	P218212

Potassium

Product Description	Grade	Form	Size	Cat. No.
Potassium (Continued)				
Potassium Citrate Monohydrate	USP/FCC	Granular	12kg	P21812
Potassium Citrate Monohydrate	USP/FCC	Granular	200 lb.	P218200LB
Potassium Citrate Monohydrate	Certified	Crystals	500g	P222500
Potassium Cyanide	Certified	Granular	100g	P2231100
Potassium Cyanide	Certified	Granular	500g	P2231500
Potassium Dichromate	Technical	Fine Granules	500g	P186500
Potassium Dichromate	Technical	Fine Granules	3kg	P1863
Potassium Dichromate	ACS Certified	Crystals	100g	P188100
Potassium Dichromate	ACS Certified	Crystals	500g	P188500
Potassium Dichromate	ACS Certified	Crystals	3kg	P1883
Potassium Ferricyanide	ACS Certified	Crystals	500g	P232500
Potassium Ferrocyanide Trihydrate	Technical	Crystals	3kg	P2343
Potassium Ferrocyanide Trihydrate	ACS Certified	Crystals	500g	P236500
Potassium Fluoride, Anhydrous	Certified	White Powder	500g	06113500
Potassium Fluoride, Anhydrous	Certified	White Powder	1kg	061131
Potassium Hydroxide	NF/FCC	Pellets	500g	P251500
Potassium Hydroxide	NF/FCC/EP/BP/JP	Pellets	2.5kg	P258212
Potassium Hydroxide	Technical	Flakes	3kg	P2463
Potassium Hydroxide	NF/FCC	Pellets	3kg	P2513
Potassium Hydroxide	NF/FCC	Pellets	50kg	P25150
Potassium Hydroxide	NF/FCC/EP/BP/JP	Pellets	50kg	P25850
Potassium Hydroxide	ACS Certified	Pellets	500g	P250500
Potassium Hydroxide	ACS Certified	Pellets	1kg	P2501
Potassium Hydroxide	ACS Certified	Pellets	3kg	P2503
Potassium Hydroxide	ACS Certified	Pellets	10kg	P25010
Potassium Hydroxide	ACS Certified	Pellets	50kg	P25050
Potassium Iodate	ACS Certified	Powder	100g	P253100
Potassium Iodate	ACS Certified	Powder	500g	P253500
Potassium Iodide	USP/FCC	Granular, Free Flowing	500g	P412500
Potassium Iodide	USP/FCC	Granular, Free Flowing	3kg	P4123
Potassium Iodide	USP/FCC	Granular, Free Flowing	10kg	P41210
Potassium Iodide	ACS Certified	Granular, Free Flowing	100g	P410100
Potassium Iodide	ACS Certified	Granular, Free Flowing	500g	P410500
Potassium Iodide	ACS Certified	Granular, Free Flowing	3kg	P4103
Potassium Iodide	ACS Certified	Granular, Free Flowing	10kg	P41010
Potassium Metabisulfite	NF/FCC	Granular	500g	P269500
Potassium Metabisulfite	NF/FCC	Granular	2.5kg	P269212
Potassium Metabisulfite	NF/FCC	Granular	12kg	P26912
Potassium Metabisulfite	Certified	Crystals	500g	P197500
Potassium Nitrate	ACS Certified	Crystals	100g	P263100
Potassium Nitrate	ACS Certified	Crystals	500g	P263500
Potassium Nitrate	ACS Certified	Crystals	3kg	P2633
Potassium Nitrate	ACS Certified	Crystals	50kg	P26350
Potassium Nitrate	Laboratory	Granular	3kg	P2613

Sodium

Product Description	Grade	Form	Size	Cat. No.
Potassium Nitrite	ACS Certified	Crystals	250g	P267250
Potassium Nitrite	ACS Certified	Crystals	500g	P267500
Potassium Oxalate Monohydrate	ACS Certified	Crystals	250g	P273250
Potassium Oxalate Monohydrate	ACS Certified	Crystals	500g	P273500
Potassium Permanganate	ACS Certified	Crystals	500g	P279500
Potassium Permanganate	ACS Certified	Crystals	2.5kg	P279212
Potassium Persulfate	ACS Certified		100g	P282100
Potassium Persulfate	ACS Certified		500g	P282500
Potassium Persulfate	Certified	Crystalline Powder	500g	P281500
Potassium Phosphate Dibasic Anhydrous	USP	Crystalline Powder	500g	P290500
Potassium Phosphate Dibasic Anhydrous	USP	Crystalline Powder	2.5kg	P290212
Potassium Phosphate Dibasic Anhydrous	USP	Crystalline Powder	12kg	P29012
Potassium Phosphate Dibasic Anhydrous	ACS Certified	Crystalline Powder	100g	P288100
Potassium Phosphate Dibasic Anhydrous	ACS Certified	Crystalline Powder	500g	P288500
Potassium Phosphate Dibasic Anhydrous	ACS Certified	Crystalline Powder	3kg	P2883
Potassium Phosphate Dibasic Anhydrous	ACS Certified	Crystalline Powder	50kg	P28850
Potassium Phosphate Dibasic Anhydrous	ACS Certified	Crystalline Powder	200 lb.	P288200LB
Potassium Phosphate Monobasic	NF/FCC	Crystals	500g	P380500
Potassium Phosphate Monobasic	NF/FCC/EP/BP	Crystals	500g	P386500
Potassium Phosphate Monobasic	NF/FCC	Crystals	2.5kg	P380212
Potassium Phosphate Monobasic	NF/FCC	Crystals	12kg	P38012
Potassium Phosphate Monobasic	NF/FCC	Crystals	250 lb.	P380250LB
Potassium Phosphate Monobasic	ACS Certified	Crystals	500g	P285500
Potassium Phosphate Monobasic	ACS Certified	Crystals	3kg	P2853
Potassium Phosphate Monobasic	ACS Certified	Crystals	10kg	P28510
Potassium Phosphate Monobasic	ACS Certified	Crystals	50kg	P28550
Potassium Phosphate Monobasic	ACS Certified	Crystals	250 lb.	P285250LB
Potassium Phosphate Monobasic	HPLC	Crystals	1kg	P2861
Potassium Phosphate Monobasic	Laboratory	Crystals/Primary Standard	500g	P382500
Potassium Phosphate Monobasic	Laboratory	Granular	500g	P284500
Potassium Phosphate Tribasic	Certified	Granular Powder	500g	P289500
Potassium Sulfate	ACS Certified	Crystals	500g	P304500
Potassium Sulfate	ACS Certified	Powder	500g	P305500
Potassium Sulfate	ACS Certified	Crystals	3kg	P3043
Potassium Sulfate	ACS Certified	Crystals	10kg	P30410
Potassium Tartrate Hemihydrate	Certified	Crystals	500g	P313500
Potassium Thiocyanate	ACS Certified	Crystals	100g	P317100
Potassium Thiocyanate	ACS Certified	Crystals	500g	P317500
Sodium				
Sodium Acetate Anhydrous	ACS Certified	Fused Crystals	500g	S210500
Sodium Acetate Anhydrous	ACS Certified	Fused Crystals	2kg	S2102
Sodium Acetate Trihydrate	USP/FCC	Granular	500g	S608500
Sodium Acetate Trihydrate	USP/FCC/EP/BP/JP	Granular	500g	S609500
Sodium Acetate Trihydrate	USP/FCC/EP/BP	White Crystals	500g	S607500
Sodium Acetate Trihydrate	USP/FCC/EP/BP/JP	Granular	2.5kg	S609212

Sodium

Product Description	Grade	Form	Size	Cat. No.
Sodium (Continued)				
Sodium Acetate Trihydrate	USP/FCC/EP/BP	White Crystals	2.5kg	S607212
Sodium Acetate Trihydrate	USP/FCC	Granular	3kg	S6083
Sodium Acetate Trihydrate	Technical	Crystals	10kg	S20710
Sodium Acetate Trihydrate	USP/FCC	Granular	12kg	S60812
Sodium Acetate Trihydrate	USP/FCC/EP/BP	White Crystals	100kg	S607100
Sodium Acetate Trihydrate	ACS Certified	Crystals	500g	S209500
Sodium Acetate Trihydrate	ACS Certified	Crystals	3kg	S2093
Sodium Acetate Trihydrate	ACS Certified	Crystals	10kg	S20910
Sodium Acetate Trihydrate	ACS Certified	Crystals	50kg	S20950
Sodium Acetate Trihydrate	HPLC		1kg	S2201
Sodium Ammonium Phosphate Tetrahydrate	Certified	Crystals	500g	S218500
Sodium Azide	Laboratory	Granular	25g	S227125
Sodium Azide	Laboratory	Granular	100g	S2271100
Sodium Azide	Laboratory	Granular	500g	S2271500
Sodium Azide	Laboratory	Granular	1kg	S22711
Sodium Benzoate	NF/FCC	Granular	500g	S224500
Sodium Bicarbonate	USP/FCC/EP/BP/JP	Crystalline Powder	2.5kg	S637212
Sodium Bicarbonate	USP/FCC	Powder	3kg	S6313
Sodium Bicarbonate	USP/FCC	Powder	10kg	S63110
Sodium Bicarbonate	USP/FCC	Powder	50kg	S63150
Sodium Bicarbonate	ACS Certified	Powder	500g	S233500
Sodium Bicarbonate	ACS Certified	Powder	3kg	S2333
Sodium Bicarbonate	ACS Certified	Powder	10kg	S23310
Sodium Bicarbonate	ACS Certified	Powder	50kg	S23350
Sodium Bicarbonate	ACS Certified	Powder	300 lb.	S233300LB
Sodium Bisulfate Monohydrate	Certified	Crystals	500g	S240500
Sodium Bisulfate Monohydrate	Certified	Crystals	3kg	S2403
Sodium Bisulfite	ACS Certified	Powder	500g	S654500
Sodium Bisulfite	ACS Certified	Powder	3kg	S6543
Sodium Bitartrate Monohydrate	Certified	Crystals	500g	S245500
Sodium Borohydride	Laboratory	Powder	10g	S67810
Sodium Borohydride	Laboratory	Powder	25g	S67825
Sodium Bromate	Certified		250g	S253250
Sodium Bromate	Certified		500g	S253500
Sodium Bromide	ACS Certified	Granular	500g	S255500
Sodium Bromide	ACS Certified	Granular	3kg	S2553
Sodium Carbonate Anhydrous	ACS Certified	Powder	500g	S263500
Sodium Carbonate Anhydrous	ACS Certified	Powder	1kg	S2631
Sodium Carbonate Anhydrous	ACS Certified	Powder	3kg	S2633
Sodium Carbonate Anhydrous	ACS Certified	Powder	10kg	S26310
Sodium Carbonate Anhydrous	ACS Certified	Powder	50kg	S26350
Sodium Carbonate Anhydrous	HPLC	Powder	500g	S495500
Sodium Carbonate Anhydrous	Laboratory	Calcined	10kg	S26110
Sodium Carbonate Monohydrate	NF/FCC/EP/BP	Granular	2.5kg	S638212

Sodium

Product Description	Grade	Form	Size	Cat. No.
Sodium Carbonate Monohydrate	NF/FCC	Crystals	10kg	S63610
Sodium Carbonate Monohydrate	ACS Certified	Crystals	3kg	S2623
Sodium Chlorate	Certified	Crystals	500g	S268500
Sodium Chloride	USP/FCC/EP/BP	Colorless to White Crystal Granules	500g	S641500
Sodium Chloride	USP/FCC/EP/BP/JP	Colorless to White Crystal Granules	500g	S642500
Sodium Chloride	USP/FCC	Granular	500g	S640500
Sodium Chloride	USP/FCC/EP/BP	Colorless to White Crystal Granules	2.5kg	S641212
Sodium Chloride	USP/FCC/EP/BP/JP	Colorless to White Crystal Granules	2.5kg	S642212
Sodium Chloride	USP/FCC	Granular	3kg	S6403
Sodium Chloride	USP/FCC	Granular	10kg	S64010
Sodium Chloride	USP/FCC/EP/BP/JP	Colorless to White Crystal Granules	12kg	S64212
Sodium Chloride	USP/FCC	Granular	50kg	S64050
Sodium Chloride	USP/FCC/EP/BP	Colorless to White Crystal Granules	350 lb.	S641P350LB
Sodium Chloride	USP/FCC/EP/BP/JP	Colorless to White Crystal Granules	350 lb.	S642350LB
Sodium Chloride	USP/FCC	Granular	350 lb.	S640350LB
Sodium Chloride	USP/FCC/EP/BP	Colorless to White Crystal Granules	1350 lb.	S641350LB
Sodium Chloride	ACS Certified	Crystals	500g	S271500
Sodium Chloride	ACS Certified	Crystals	1kg	S2711
Sodium Chloride	ACS Certified	Crystals	3kg	S2713
Sodium Chloride	ACS Certified	Crystals	10kg	S27110
Sodium Chloride	ACS Certified	Crystals	50 lb.	S271350LB
Sodium Chloride	ACS Certified	Crystals	50kg	S27150
Sodium Chloride	Certified	Crystals/Biological	500g	S671500
Sodium Chloride	Certified	Crystals/Biological	3kg	S6713
Sodium Chloride	Certified	Crystals/Biological	10kg	S67110
Sodium Chromate Anhydrous	Certified	Powder	250g	S272250
Sodium Chromate Anhydrous	Certified	Powder	500g	S272500
Sodium Citrate Dihydrate	USP/FCC/EP/BP/JP	Granular	2.5kg	S470212
Sodium Citrate Dihydrate	USP/FCC	Granular	3kg	S4663
Sodium Citrate Dihydrate	USP/FCC	Powder	3kg	S4673
Sodium Citrate Dihydrate	Certified	Granular	500g	S279500
Sodium Citrate Dihydrate	Certified	Granular	3kg	S2793
Sodium Citrate Dihydrate	Certified	Granular	10kg	S27910
Sodium Citrate Dihydrate	Certified	Granular	50kg	S27950
Sodium Citrate Dihydrate	Certified	Granular	275 lb.	S279275LB
Sodium Cobaltinitrite	ACS Certified	Powder	100g	S280100
Sodium Cyanide	Technical	Granular	250g	S283I250
Sodium Cyanide	Technical	Granular	500g	S283I500
Sodium Cyanide	ACS Certified	Granular	100g	S284I100
Sodium Cyanide	ACS Certified	Granular	500g	S284I500
Sodium Dichromate	ACS Certified	Crystals	3kg	S2583
Sodium Dichromate Dihydrate	Technical	Granular	500g	S234500
Sodium Dichromate Dihydrate	Technical	Granular	3kg	S2343
Sodium Dichromate Dihydrate	Technical	Granular	10kg	S23410
Sodium Dichromate Dihydrate	ACS Certified	Crystals	100g	S235100

Sodium

Product Description	Grade	Form	Size	Cat. No.
Sodium (Continued)				
Sodium Dichromate Dihydrate	ACS Certified	Crystals	500g	S235500
Sodium Dichromate Dihydrate	ACS Certified	Crystals	3kg	S2353
Sodium Diethyldithiocarbamate Trihydrate	ACS Certified	Crystals	100g	S287100
Sodium Dodecyl Sulfate	Certified		25g	0267425
Sodium Fluoride	ACS Certified	Powder	100g	S299100
Sodium Fluoride	ACS Certified	Powder	500g	S299500
Sodium Fluoride	ACS Certified	Powder	3kg	S2993
Sodium Formate	ACS Certified	Crystals	500g	S648500
Sodium Hydrosulfite	Laboratory	Powder	100g	S310100
Sodium Hydrosulfite	Laboratory	Powder	500g	S310500
Sodium Hydroxide	NF/EP/BP/FCC		2.5kg	S399212
Sodium Hydroxide	NF/EP/BP/FCC		500g	S399500
Sodium Hydroxide	NF/FCC	Pellets	500g	S320500
Sodium Hydroxide	NF/EP/BP/FCC		1kg	S3991
Sodium Hydroxide	NF/FCC	Pellets	1kg	S3201
Sodium Hydroxide	NF/FCC/EP/BP/JP	Granular	2.5kg	S392212
Sodium Hydroxide	NF/FCC	Pellets	3kg	S3203
Sodium Hydroxide	NF/FCC	Pellets	10kg	S32010
Sodium Hydroxide	NF/FCC/EP/BP/JP	Granular	12kg	S39212
Sodium Hydroxide	NF/EP/BP/FCC		50kg	S39950
Sodium Hydroxide	NF/FCC/EP/BP/JP	Granular	50kg	S39250
Sodium Hydroxide	NF/FCC	Pellets	50kg	S32050
Sodium Hydroxide	ACS Certified	Pellets	100g	S318100
Sodium Hydroxide	ACS Certified	Pellets	500g	S318500
Sodium Hydroxide	ACS Certified	Pellets	1kg	S3181
Sodium Hydroxide	ACS Certified	Pellets	3kg	S3183
Sodium Hydroxide	ACS Certified	Pellets	5kg	S3185
Sodium Hydroxide	ACS Certified	Pellets	10kg	S31810
Sodium Hydroxide	ACS Certified	Pellets	50kg	S31850
Sodium Hydroxide	Laboratory	Beads	3kg	S6123
Sodium Hydroxide	Laboratory	Beads	3kg	S6133
Sodium Hydroxide	Laboratory	Beads	10kg	S61310
Sodium Hydroxide	Laboratory	Beads	50kg	S61250
Sodium Hydroxide	Laboratory	Beads	50kg	S61350
Sodium Iodate	Certified	Powder	100g	S322100
Sodium Iodide	USP	Granular	500g	S473500
Sodium Iodide	Certified	Crystals	100g	S324100
Sodium Iodide	Certified	Crystals	500g	S324500
Sodium L-(+)-Tartrate Dihydrate	Certified	Crystals	100g	S433100
Sodium L-(+)-Tartrate Dihydrate	Certified	Crystals	500g	S433500
Sodium Lauryl Sulfate	NF/FCC	Powder	500g	S529500
Sodium Lauryl Sulfate	NF/FCC	Powder	3kg	S5293
Sodium Metabisulfite	NF/FCC	Granular	500g	S242500
Sodium Metabisulfite	NF/FCC	Granular	2.5kg	S242212

Sodium

Product Description	Grade	Form	Size	Cat. No.
Sodium Metabisulfite	NF/FCC	Granular	12kg	S24212
Sodium Metabisulfite	ACS Certified	Powder	500g	S244500
Sodium Metabisulfite	ACS Certified	Powder	3kg	S2443
Sodium Metabisulfite	Laboratory		10kg	S24310
Sodium Metaperiodate	ACS Certified	Crystals	50g	S39850
Sodium Metaperiodate	ACS Certified	Crystals	100g	S398100
Sodium Metaphosphate	Laboratory	Powder	500g	S333500
Sodium Metasilicate Nonahydrate	Certified	Crystals	500g	S408500
Sodium Metasilicate Nonahydrate	Certified	Crystals	3kg	S4083
Sodium Metasilicate Pentahydrate	Technical	Crystals	3kg	S4073
Sodium Methoxide	Laboratory	Powder	100g	S335100
Sodium Molybdate Dihydrate	ACS Certified	Crystalline Powder	500g	S336500
Sodium Molybdate Dihydrate	ACS Certified	Crystalline Powder	3kg	S3363
Sodium Nitrate	ACS Certified	Crystals	500g	S343500
Sodium Nitrate	ACS Certified	Crystals	3kg	S3433
Sodium Nitrate	Laboratory	Granular	3kg	S3423
Sodium Nitrite	USP/FCC	Flakes	3kg	S3383
Sodium Nitrite	ACS Certified	Crystals	250g	S347250
Sodium Nitrite	ACS Certified	Crystals	500g	S347500
Sodium Nitrite	ACS Certified	Crystals	3kg	S3473
Sodium Nitrite	ACS Certified	Crystals	10kg	S34710
Sodium Nitroferricyanide Dihydrate	ACS Certified	Crystals	100g	S350100
Sodium Oxalate	ACS Certified		500g	S487500
Sodium Perchlorate Monohydrate	HPLC		100g	S490100
Sodium Perchlorate Monohydrate	HPLC		500g	S490500
Sodium Perchlorate Monohydrate	Laboratory	Crystals	500g	S360500
Sodium Perchlorate Monohydrate	Laboratory	Crystals	2.5kg	S360212
Sodium Peroxide	ACS Certified		500g	S363500
Sodium Persulfate	Certified	White Crystals, Solid	500g	06114500
Sodium Persulfate	Certified	White Crystals, Solid	1kg	061141
Sodium Phosphate Dibasic Anhydrous	USP		500g	S375500
Sodium Phosphate Dibasic Anhydrous	USP		2.5kg	S375212
Sodium Phosphate Dibasic Anhydrous	USP		12kg	S37512
Sodium Phosphate Dibasic Anhydrous	USP/FCC/EP		12kg	S37912
Sodium Phosphate Dibasic Anhydrous	USP		250 lb.	S375250LB
Sodium Phosphate Dibasic Anhydrous	USP/FCC/EP		2.5kg	S379212
Sodium Phosphate Dibasic Anhydrous	ACS Certified	Granular or Powder	500g	S374500
Sodium Phosphate Dibasic Anhydrous	ACS Certified	Granular or Powder	1kg	S3741
Sodium Phosphate Dibasic Anhydrous	ACS Certified	Granular or Powder	3kg	S3743
Sodium Phosphate Dibasic Anhydrous	ACS Certified	Granular or Powder	50kg	S37450
Sodium Phosphate Dibasic Anhydrous	ACS Certified	Granular or Powder	250 lb.	S374250LB
Sodium Phosphate Dibasic Anhydrous	Laboratory	Powder	3kg	S3933
Sodium Phosphate Dibasic Dihydrate	USP/FCC/EP/BP		500g	S472500
Sodium Phosphate Dibasic Dihydrate	USP/FCC/EP/BP		2.5kg	S472212
Sodium Phosphate Dibasic Heptahydrate	USP	Granular	3kg	S4713

Sodium

Product Description	Grade	Form	Size	Cat. No.
Sodium (Continued)				
Sodium Phosphate Dibasic Heptahydrate	ACS Certified	Crystals	500g	S373500
Sodium Phosphate Dibasic Heptahydrate	ACS Certified	Crystals	3kg	S3733
Sodium Phosphate Dibasic Heptahydrate	ACS Certified	Crystals	250 lb.	S373250LB
Sodium Phosphate Monobasic Anhydrous	USP		500g	S397500
Sodium Phosphate Monobasic Dihydrate	Certified	Crystals	500g	S381500
Sodium Phosphate Monobasic Dihydrate	Certified	Crystals	3kg	S3813
Sodium Phosphate Monobasic Monohydrate	USP/FCC	Crystals or Granular	500g	S468500
Sodium Phosphate Monobasic Monohydrate	ACS Certified	Crystals	50g	S36950
Sodium Phosphate Monobasic Monohydrate	ACS Certified	Crystals	500g	S369500
Sodium Phosphate Monobasic Monohydrate	ACS Certified	Crystals	1kg	S3691
Sodium Phosphate Monobasic Monohydrate	ACS Certified	Crystals	3kg	S3693
Sodium Phosphate Tribasic Dodecahydrate	Technical	Granular	3kg	S3763
Sodium Phosphate Tribasic Dodecahydrate	Technical	Granular	50kg	S37650
Sodium Phosphate Tribasic Dodecahydrate	ACS Certified	Crystals	500g	S377500
Sodium Phosphate Tribasic Dodecahydrate	ACS Certified	Crystals	3kg	S3773
Sodium Phosphate Tribasic Dodecahydrate	ACS Certified	Crystals	50kg	S37750
Sodium Potassium Tartrate	USP	Powder	500g	S386500
Sodium Potassium Tartrate	USP	Powder	2.5kg	S386212
Sodium Potassium Tartrate	USP	Powder	12kg	S38612
Sodium Potassium Tartrate Tetrahydrate	ACS Certified	Crystals	500g	S387500
Sodium Potassium Tartrate Tetrahydrate	ACS Certified	Crystals	3kg	S3873
Sodium Potassium Tartrate Tetrahydrate	ACS Certified	Crystals	10kg	S38710
Sodium Pyrophosphate Decahydrate	ACS Certified	Crystals	500g	S390500
Sodium Pyrophosphate Decahydrate	ACS Certified	Crystals	3kg	S3903
Sodium Pyrophosphate Decahydrate	ACS Certified	Crystals	10kg	S39010
Sodium Salicylate	USP	Small Crystals, Free Flowing	12kg	S39612
Sodium Succinate Hexahydrate	Laboratory	Crystals	500g	S413500
Sodium Sulfate Anhydrous	USP/FCC	Powder	500g	S429500
Sodium Sulfate Anhydrous	USP/FCC	Powder	2.5kg	S429212
Sodium Sulfate Anhydrous	Technical	Low Nitrogen	3kg	S4203
Sodium Sulfate Anhydrous	Technical	Low Nitrogen	10kg	S42010
Sodium Sulfate Anhydrous	USP/FCC	Powder	12kg	S42912
Sodium Sulfate Anhydrous	USP/FCC/EP/BP	Powder	12kg	S43412
Sodium Sulfate Anhydrous	USP/FCC	Powder	250 lb.	S429250LB
Sodium Sulfate Anhydrous	ACS Certified	Granular	500g	S421500
Sodium Sulfate Anhydrous	ACS Certified	Granular, 10-60 Mesh	500g	S415500
Sodium Sulfate Anhydrous	ACS Certified	Granular	1kg	S4211
Sodium Sulfate Anhydrous	ACS Certified	Granular, 10-60 Mesh	1kg	S4151
Sodium Sulfate Anhydrous	ACS Certified	Granular, 10-60 Mesh	2.5kg	S415212
Sodium Sulfate Anhydrous	ACS Certified	Granular	3kg	S4213
Sodium Sulfate Anhydrous	ACS Certified	Granular	10kg	S42110
Sodium Sulfate Anhydrous	ACS Certified	Granular, 10-60 Mesh	10kg	S41510S
Sodium Sulfate Anhydrous	ACS Certified	Granular	50kg	S42150
Sodium Sulfate Anhydrous	ACS Certified	Granular, 10-60 Mesh	200 lb.	S415200LB

Strontium

Product Description	Grade	Form	Size	Cat. No.
Sodium Sulfate Anhydrous	ACS Certified	Granular	300 lb.	S421300LB
Sodium Sulfate Decahydrate	USP	Granular	500g	S414500
Sodium Sulfate Decahydrate	USP	Granular	2.5kg	S414212
Sodium Sulfate Decahydrate	USP	Granular	12kg	S41412
Sodium Sulfate Decahydrate	Certified	Crystals	500g	S419500
Sodium Sulfhydrate	Laboratory	Flakes	500g	S423500
Sodium Sulfide	Technical	Flakes	2.5kg	S426212
Sodium Sulfide Nonahydrate	Certified	Crystals	500g	S425500
Sodium Sulfide Nonahydrate	Certified	Crystals	2.5kg	S425212
Sodium Sulfite Anhydrous	ACS Certified	Crystals	500g	S430500
Sodium Sulfite Anhydrous	ACS Certified	Crystals	3kg	S4303
Sodium Sulfite Anhydrous	ACS Certified	Crystals	10kg	S43010
Sodium Sulfite Anhydrous	Laboratory	Crystals/Low Phosphate	500g	S447500
Sodium Sulfite Anhydrous	Laboratory	Crystals/Low Phosphate	3kg	S4473
Sodium Tetraborate Anhydrous	Laboratory	Fused, Ground	10kg	S25210
Sodium Tetraborate Decahydrate	NF	Powder or Granular	500g	S246500
Sodium Tetraborate Decahydrate	NF	Powder or Granular	2.5kg	S246212
Sodium Tetraborate Decahydrate	NF	Powder or Granular	12kg	S24612
Sodium Tetraborate Decahydrate	ACS Certified	Crystals	500g	S248500
Sodium Tetraborate Decahydrate	ACS Certified	Powder	500g	S249500
Sodium Tetraborate Decahydrate	ACS Certified	Crystals	3kg	S2483
Sodium Tetraborate Decahydrate	ACS Certified	Crystals	10kg	S24810
Sodium Tetraphenyl Boron	ACS Certified	Crystals	25g	S65225
Sodium Thiocyanate	ACS Certified	Crystals	500g	S441500
Sodium Thiosulfate	USP/FCC	Crystals	500g	S474500
Sodium Thiosulfate	USP/FCC	Crystals	3kg	S4743
Sodium Thiosulfate	USP/FCC	Crystals	12kg	S47412
Sodium Thiosulfate Anhydrous	Certified		500g	S446500
Sodium Thiosulfate Anhydrous	Certified		3kg	S4463
Sodium Thiosulfate Pentahydrate	USP/FCC/EP/BP	Crystals	2.5kg	S475212
Sodium Thiosulfate Pentahydrate	USP/FCC/EP/BP	Crystals	50kg	S47550kg
Sodium Thiosulfate Pentahydrate	USP/FCC/EP/BP	Crystals	12kg	S47512
Sodium Thiosulfate Pentahydrate	ACS Certified	Crystals	500g	S445500
Sodium Thiosulfate Pentahydrate	ACS Certified	Crystals	3kg	S4453
Sodium Thiosulfate Pentahydrate	ACS Certified	Crystals	10kg	S44510
Sodium Thiosulfate Pentahydrate	ACS Certified	Crystals	50kg	S44550
Sodium Triphosphate	Laboratory	Granular	500g	S645500
Sodium Tungstate Dihydrate	ACS Certified	Crystals	500g	S450500
Sodium Tungstate Dihydrate	ACS Certified	Crystals	3kg	S4503
Sodium Vanadate	Laboratory		50g	S45450
Strontium				
Strontium Chloride Hexahydrate	ACS Certified	Crystals	500g	S541500
Strontium Nitrate Anhydrous	ACS Certified	Crystals	100g	S549100
Strontium Nitrate Anhydrous	ACS Certified	Crystals	500g	S549500

To place an order, contact your local Fisher Scientific office.



©2014 Thermo Fisher Scientific Inc. All rights reserved.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

AMERICAS

Canada
Fisher Scientific Canada
www.fishersci.ca
Tel: 800-234-7437
Fax: 800-463-2996

Latin America
Fisher Scientific Global Export,
Latin America
www.fishersci.com
Tel: 770-871-4725
Fax: 770-871-4726

United States
Fisher HealthCare
www.fisherhealthcare.com
Tel: 800-640-0640
Fax: 800-290-0290

Fisher Scientific
www.fishersci.com
Tel: 800-766-7000
Fax: 800-926-1166

EUROPE

Austria
info.austria@thermofisher.com
www.at.fishersci.com
Tel: 0800 20 88 40
Fax: 0800 20 66 90

Belgium
be.fisher@thermofisher.com
www.be.fishersci.com
Tel: 056 260 260
Fax: 056 260 270

Czech Republic
info.cz@thermofisher.com
www.thermofisher.cz
Tel: 466 798 230
Fax: 466 435 008

Denmark
kundeservice@thermofisher.com
www.dk.fishersci.com
Tel: 70 27 99 20
Fax: 70 27 99 29

Finland
fisher.fi@thermofisher.com
www.fi.fishersci.com
Tel: 09-802 76 280
Fax: 09-802 76 235

France
fr.fisher@thermofisher.com
www.fr.fishersci.com
Tel: 03 88 67 14 14
Fax: 03 88 67 11 68

Germany
info.germany@thermofisher.com
www.de.fishersci.com
Tel: 2304 932-5
Fax: 2304 932-950

Ireland
fsie.sales@thermofisher.com
www.ie.fishersci.com
Tel: 01 885 5854
Fax: 01 899 1855

Italy
it.fisher@thermofisher.com
www.it.fishersci.com
Tel: 02 950 59 478
Fax: 02 950 59 479

The Netherlands
nl.info@thermofisher.com
www.nl.fishersci.com
Tel: 020 487 70 00
Fax: 020 487 70 70

Norway
fisher.no@thermofisher.com
www.no.fishersci.com
Tel: 22 95 59 59
Fax: 22 95 59 40

Portugal
pt.fisher@thermofisher.com
www.pt.fishersci.com
Tel: 21 425 33 50
Fax: 21 425 33 51

Spain
es.fisher@thermofisher.com
www.es.fishersci.com
Tel: 902 239 303
Fax: 902 239 404

Sweden
fisher.se@thermofisher.com
www.se.fishersci.com
Tel: 31 352 32 00
Fax: 31 352 32 50

Switzerland
info.ch@thermofisher.com
www.ch.fishersci.com
Tel: 056 618 41 11
Fax: 056 618 41 41

United Kingdom
fsuk.sales@thermofisher.com
www.fisher.co.uk
Tel: 01509 555500
Fax: 01509 555111

ASIA

India
qfc.customer@thermofisher.com
www.fishersci.in
Tel: +91-22-6680 3000
Fax: +91-22-6680 3001

Korea
sales.fsk@thermofisher.com
www.fishersci.co.kr
Tel: 02-3420-8700
Fax: 02-3420-8710

Malaysia
enquiry.my@thermofisher.com
www.fishersci.com.my
Tel: (603) 51228888
Fax: (603) 5121 8899

Singapore
enquiry.sg@thermofisher.com
www.fishersci.com.sg
Tel: (65) 6873 6006
Fax: (65) 6873 5005

OCEANIA

Australia
AUinfo@thermofisher.com
www.thermofisher.com.au
Tel: 1300-735-292
Fax: 1800-067-639

New Zealand
NZinfo@thermofisher.com
www.thermofisher.co.nz
Tel: 0800-933-966
Fax: +64 9 980 6788

