

# WHEATON

Buy 2 Cases of BOD Bottles, and get a FREE BOD Bottle Rack\*



## BOD Bottles

- Ideal for incubating diluted samples of sewage, sewage effluents, polluted waters, and industrial wastes to determine the amount of oxygen required during the stabilization of the decomposable organic matter by aerobic biochemical action
- We recommend the 300mL capacity bottles for the five-day BOD test as referenced in Standard Methods for the Examination of Water and Wastewater, 21st Edition, 5210B, 2005
- Specially designed bottle shoulder radius that sweeps all air from inside the bottle during filling
- The interchangeable stoppers have a tapered bottom that prevents air entrapment
- The bottles have a flared mouth to form a water seal around the stopper that prevents air from being drawn into the bottle during incubation
- Stopper joint is compatible with the probes of the leading meters for BOD and dissolved oxygen
- Large, permanently screened-on writing patch on all bottles
- Manufactured from USP Type I borosilicate glass

Cat. No.	Number Sequence	Barcoding	Dia x H (mm)	Stopper	Qty/Case
<b>60mL Capacity</b>					
02926100A	Un-numbered	No	43 x 115	Glass Pennyhead	36
02-926-12	Un-numbered	No	43 x 115	Glass Robotic	36
02926100B	01-36	Yes	43 x 115	Glass Pennyhead	36
02-926-13	01-36	Yes	43 x 115	Glass Robotic	36
02926100C	37-72	Yes	43 x 115	Glass Pennyhead	36
02-926-14	37-72	Yes	43 x 115	Glass Robotic	36
<b>300mL Capacity</b>					
02-926-7A	Un-numbered	No	69 x 165	Glass Pennyhead	24
02-926-26	Un-numbered	No	69 x 165	Glass Robotic	24
02-926-8	01-24	Yes	69 x 165	Glass Pennyhead	24
02-926-27	01-24	Yes	69 x 165	Glass Robotic	24
02-926-9	25-48	Yes	69 x 165	Glass Pennyhead	24
02-926-28	25-48	Yes	69 x 165	Glass Robotic	24
02-926-11	49-72	Yes	69 x 165	Glass Pennyhead	24
02-926-29	49-72	Yes	69 x 165	Glass Robotic	24
02-926-25	73-96	Yes	69 x 165	Glass Pennyhead	24
02-926-30	73-96	Yes	69 x 165	Glass Robotic	24
02-926-90	97-120	Yes	69 x 165	Glass Pennyhead	24
02-926-31	97-120	Yes	69 x 165	Glass Robotic	24

\* For additional number sequences or specials contact WHEATON Customer Service.

## BOD Bottle Rack

- Conveniently and safely transports BOD bottles
- Racks are stackable

Cat. No.	For Bottle Size	Holds	L x W x H (in)*	Qty/Case
02-927-031	60mL	20 Bottles	13 x 10 x 7 in	1
02-927-030	300mL	12 Bottles	13 x 10 x 7 in	1

\*33.03 x 25.4 x 17.78cm

Complete the redemption coupon below and send or mail/fax/email it along with your proof of purchase to:

**Limited Time Offer**  
Offer Valid thru December 31, 2012

### WHEATON

Attention: Joel Torres  
1501 North 10th Street  
Millville, NJ 08332-2038 USA

Email: joel.torres@wheaton.com  
800.225.1437 (T) | 856.825.9035 (F)

Please send my FREE BOD Bottle Rack,  
Cat. No. \_\_\_\_\_.

Please enclose or fax a copy of your proof of purchase: Invoice, packing slip, or purchase order with your redemption.

Name: \_\_\_\_\_  
Company: \_\_\_\_\_  
Address: \_\_\_\_\_  
Address: \_\_\_\_\_  
City: \_\_\_\_\_  
State: \_\_\_\_\_  
Zip/Postal: \_\_\_\_\_  
Country: \_\_\_\_\_  
Tel: \_\_\_\_\_  
E-Mail: \_\_\_\_\_

This offer is void where prohibited by law or company policy. All information must be completed to be a valid entry. Offer valid in North America only thru 12/31/2012. Limited to two redemptions per customer. \*Not for resale.