

BD Falcon™ 50 ml Conical Tube with Flip-Top Cap

New innovative cap design—*same tube you have always trusted.*



The BD Falcon™ 50 ml Conical Tube with Flip-Top Cap is ideal for any research applications requiring one-handed operations to open and close the cap. This novel design saves time and effort in applications requiring multiple aliquoting, storage, and pouring from the same tube, while maintaining the same superior quality and performance of our standard screw-cap closure.

This new addition to our BD Falcon Conical Tube product family offers a number of unique features and benefits:

- **Easy-to-Use**—One-hand cap manipulation is virtually effortless with the unique hinge action of the flip-top cap. Easy-to-read blue graduations on the tube allow rapid volumetric assessment.
- **Splash-Guard**—Innovative design of the cap allows opening and closing of the cap without risk of splatter.
- **Secure Seal**—Innovative click-bead ensures a tight seal every time.
- **Ambidextrous**—Thumb grips on each side of the cap ensure a sterile opening technique.
- **High Strength**—The BD Falcon Conical Tube with Flip-Top Cap can be centrifuged up to 9400 x g.
- **Non-Toxic**—Resins are selected via U.S. Pharmacopeia (USP) toxicity tests to ensure our tubes do not contain concentrations of bioreactive substances that could potentially interfere with biological experiments.
- **Optical Clarity**—High-clarity polypropylene conical tubes provide easy viewing through ultraclear tube walls.

BD Falcon™ 50 ml Conical Tube with Flip-Top Cap



BD Falcon™ 50 ml Conical Tubes with Flip-Top Cap offer:

- Easy-to-open with one hand
- Secure and leak resistant cap
- Minimizes accidental cap mix-up
- Sterile, non-pyrogenic
- Fits standard centrifuge rotors
- R.C.F. rating of 9400 x g
- Large white, solvent-resistant writing patch
- Reduces unnecessary environmental exposure
- Class I Exempt Medical Device (US), General Laboratory Use (OUS)
- Medical-grade packaging



Ideally suited for serial dilutions and the safe storage of solutions, the BD Falcon 50 ml Conical Tube with Flip-Top Cap is a valuable addition to the BD Falcon product family. Visit bdbiosciences.com/fliptop to learn more about how this product will bring added convenience to your life science research.

Ordering Information

Fisher Cat. No.	BD Cat. No.	Description	Qty./Pack	Qty./Case
BD Falcon™ 50 ml Conical Tube with Flip-Top Cap				
02-683-173	352077	BD Falcon™ 50 ml Conical Tube with Flip-Top Cap	22	440

To place an order in the U.S., contact Fisher Scientific at **800.766.7000**.

For technical assistance, contact Technical Support at:
tel: **877.232.8995** or 978.901.7389; fax: 978.901.7491; e-mail: labware@bd.com



**Fisher
Scientific**

For customer service, call **1-800-766-7000**.
To fax an order, use **1-800-926-1166**.
To order online: www.fishersci.com