

# Fisher Scientific Rockers and Shakers

Free accessory with purchase

**Purchase a Fisher Scientific Microplate Shaker and receive a 4 place plate top FREE when you use Promo No. 13-687-708PM**

## Fisher Scientific™ Digital Microplate Shaker

Compact, variable speed

**These compact, horizontal shakers deliver reliable and regulated shaking for up to four microplates.**

- Variable shaking speed: 150 to 1,000rpm
- Quick and easy fitting of any standard-depth multiwall plates
- Continuous or timed operation, with automatic switch off
- Operating temperature 4° to 40°C allows for use in cold rooms or (non CO<sub>2</sub>) incubators
- Easy-to-use digital timer ensures accurate count-down and repeatability of time sensitive incubations

### Free accessory for Fisher Scientific™ Digital Microplate Shaker

- Cat. No: 13-687-709
- 4-plate platform
- 6.6 x 10.3" platform

**Purchase a Fisher Scientific Digital Multi-Platform Shaker and receive a FREE flat adjustable bar platform with rubber non-slip mat when you use Promo No. 13-687-700PM**

## Fisher Scientific™ Digital Multi-Platform Shaker

Easy-to-use, microprocessor-controlled

**This multi-platform shaker provides a smooth, quiet orbital motion for mixing in bottles, flasks and beakers.**

- Variable shaking speed: 50 to 450rpm
- Holds up to 8.8lbs
- Automatic load balancing system
- Continuous or timed operation with automatic switch-off
- Maximum flexibility with an option of four interchangeable platforms
- Simple set up and easy to operate
- Low voltage power supply provides for safe, cold room operation and low-energy consumption
- Operating temperature range of 4° to 40°C allows for use in cold rooms or (non CO<sub>2</sub>) incubators

### Universal Platform for for Fisher Scientific™ Digital Multi-Platform Shaker


(SHOWN ON UNIT)

- Cat No: 13-687-710
- Dimensions: 13.4 x 16.8 x 4.1 in
- Universal platform with adjustable bars



13-687-708

Use Promo No. **13-687-708PM** to receive free 4 top plate (pictured below)



13-687-709



Unit Cat No.:  
**13-687-700**  
Platform Cat. No.:  
**13-687-710**

Use Promo No. **13-687-700PM** to receive unit with free platform



Designed and tested to comply with U.S., Canada and European safety standards

# Fisher Scientific Rockers and Shakers

Free accessory with purchase

**Purchase a Fisher Scientific Digital 3D Platform Rotator, Digital Multi-Function 3D Rotator, or Platform Rocker and receive a FREE dimpled mat (shown below) when you use the promotional part number when ordering.**

## Fisher Scientific™ Rockers and Rotators

Particularly suited for routine laboratory applications

**These compact rockers and rotators provide smooth action for gentle sample agitation.**

- Ideal for tubes, culture flasks, dishes and boxes
- Fits neatly into the smallest workspace
- Standard unit comes with flat top with no slip mat
- Operating temperature range of 4° to 40°C allows for use in cold rooms or (non CO<sub>2</sub>) incubators
- Each rocker/rotator has unique features and benefits, please see website or product brochure for full listing of these additional features

## Dimpled mat for Fisher Scientific Rockers and Rotators

- Ideal for mixing samples in any size tube
- 7.8 x 7.8" dimpled mat



13-687-706



**Fisher Scientific  
Platform Rocker  
13-687-704**

Use Promo No. **13-687-704PM** to receive free dimpled mat with unit



**Fisher Scientific  
Multi-Function 3D Rotator  
13-687-707**

Use Promo No. **13-687-707PM** to receive free dimpled mat with unit



**Fisher Scientific  
3D Platform Rotator  
13-687-705**

Use Promo No. **13-687-705PM** to receive free dimpled mat with unit

Must be purchased using the promotional catalog number listed to take advantage of special offers. Offers expire December 31, 2014.

© 2014 Thermo Fisher Scientific Inc. All rights reserved.  
Trademarks used are owned as indicated at [www.fishersci.com/trademarks](http://www.fishersci.com/trademarks).

### In the United States:

For customer service, call 1-800-766-7000  
To fax an order, use 1-800-926-1166  
To order online: [www.fishersci.com](http://www.fishersci.com)

13\_2015 JA/JC 12/13 BN1210131

### In Canada:

For customer service, call 1-800-234-7437  
To fax an order, use 1-800-463-2996  
To order online: [www.fishersci.ca](http://www.fishersci.ca)



A Thermo Fisher Scientific Brand